



GeneAll®



2012 / 2013 GeneAll
Innovative Life Science System

GeneAll[®] Innovative Life Science System

With the advance in molecular biological techniques, researchers have preferred the commercial ready-made kits to lab-made reagents in order to concentrate on doing research itself rather than making reagents.

GeneAll[®] DNA and RNA Purification kit series are basic materials in molecular biological experiments and offer fast, accurate, convenient and reproducible methods. Every GeneAll[®] product is manufactured under strictly clean condition and controlled thoroughly from lot to lot, and we proudly guarantee the stable and consistent quality.

GeneAll[®] SV column contains silica membrane that will bind DNA and easily apply to both centrifugation and vacuum protocols. Purification step is so simple, bind-wash-elute, that is all. Under high salt condition, DNA bind to silica membrane and impurities pass through membrane into a collection tube. The membranes are washed with an ethanol-containing buffer to remove any residual of proteins, cellular debris, salts, remnant of agarose, enzymatic reaction components and etc. Finally DNA is released into a clean collection tube with water or low ionic strength buffer.

GeneAll[®] Plasmid DNA Purification System



GeneAll[®] Plasmid DNA Purification Systems utilize glass microfiber membranes based on the modified alkaline lysis method. GeneAll[®] Hybrid-Q[™] Plasmid Rapidprep with new patented EzClear[™] filter provides the alternative methods for standard or rapid preparation of Plasmid DNA depending on plasmid copy number, host strain, culture medium and culture volume.

GeneAll[®] Exfection[™] Plasmid LE (Low Endotoxin) and EF (Endotoxin-Free) kits provide simple and fast method for the purification of

plasmid DNA with low endotoxin contaminants. Endotoxins present in the cell membrane of gram-negative bacteria are common contaminants in plasmid preparations and can significantly reduce transfection efficiencies. Exfection[™] series can be used for the transfection of most cell lines through the removal of endotoxins: advanced phase separation and endotoxin removal washing.

GeneAll[®] Fragment DNA Purification System



GeneAll[®] Expin[™] series provide reliable and fast methods for the purification of fragment DNA from agarose gel and PCR or enzymatic reaction mixtures. Expin[™] Gel SV takes advantage of glass fiber membrane to recover DNA of 80 bp to 10 kb from most grades of agarose gel in yields reaching 90 %. Expin[™] PCR SV is used to recover DNA of 100 bp to 10 kb from PCR or enzymatic reaction mixtures and very effective to the removal of PCR primer dimer. Expin[™] Combo GP kit is the combined product of Expin[™] Gel SV and Expin[™]

PCR SV. Expin[™] CleanUp SV is designed for fast and simple method for purification of fragment DNA of 40 bp to 10 kb from various enzymatic reactions in just 5 minutes.



GeneAll® Genomic DNA Purification System



GeneAll® Exgene™ and **GenEx™** series are designed for the purification of total DNA from a variety of sample sources. Exgene™ series provide fast and easy method in convenient spin or vacuum column format and are no need phenol extraction or alcohol precipitation.

GeneAll® GenEx™ series provide convenient, scalable purification method in the specially formulated buffer systems. Purified total DNA can be directly applicable in conventional PCR, real-time PCR, Southern blotting, genotyping, RFLP and other downstream applications.

GeneAll® DirEx™ can be conveniently used for preparation of the total DNA from various biological samples without the use of toxic chemical such as phenol or chloroform.

GeneAll® RNA Purification System

GeneAll® RiboEx™ series are designed for total RNA isolation from various samples. RiboEx™ is based on the disruption of cells in a monophasic lysis solution containing phenol and salt followed by alcohol precipitation of the RNA. Hybrid-R™ eliminates alcohol precipitation by binding of RNA with column. RiboEx™ LS is a concentrated form of RiboEx™ for total RNA isolation from liquid samples, while RiboEx™ is more suitable for solid samples and pelleted cells. Riboclear™ provides an easy and rapid method for

RNA cleanup or concentration from various RNA samples in just 6 minutes. Ribospin™ series provide fast and easy method in convenient spin column format and isolate highly purified RNA in 15 minutes. Allspin™ provides the method for the isolation of genomic DNA and total RNA simultaneously from a single sample of tissue or cultured cells.



GeneAll® PCR Amplification System



GeneAll® AmpONE™ PCR enzyme series are designed for various applications such as cloning, genotyping, sequencing, routine PCR, real-time PCR and long PCR. It is provided with HQ buffer for the amplification of a higher order structure. AmpONE™ Premix (lyophilized or solution form) and AmpMaster™ mix (2-fold liquid form) are provided with DNA polymerase, reaction buffer, dNTP and loading dye.

GeneAll[®] Table of Contents

Plasmid DNA Purification System

Page

Hybrid-Q [™] Plasmid Rapidprep	08
Exprep [™] Plasmid SV	10
Exfection [™] Plasmid LE/EF	12
Selection Guide for Plasmid DNA Purification	14

Fragment DNA Purification System

Expin [™] Gel SV	18
Expin [™] PCR SV	19
Expin [™] CleanUp SV	20
Expin [™] Combo GP	21
Selection Guide for Fragment DNA Purification	22

Genomic DNA Purification System

Exgene [™] Tissue SV (<i>plus!</i>)	24
Exgene [™] Blood SV	26
Exgene [™] Clinic SV	28
Exgene [™] Cell SV	30
Exgene [™] Plant SV	32
Exgene [™] Soil DNA mini	34
Exgene [™] Genomic DNA micro	36
Exgene [™] GMO SV	38
GenEx [™] Blood / Cell / Tissue	40
GenEx [™] Plant (<i>plus!</i>)	42
DirEx [™]	44
Selection Guide for Genomic DNA Purification	46

GeneAll[®] Table of Contents

RNA Purification System

Page

Hybrid-R [™]	50
Hybrid-R [™] Blood RNA	52
Hybrid-R [™] miRNA	54
RiboEx [™]	56
RiboEx [™] LS	58
Ribospin [™]	60
Ribospin [™] vRD (<i>plus!</i>)	62
Ribospin [™] Plant	64
Riboclear [™] (<i>plus!</i>)	66
Allspin [™]	68
Selection Guide for RNA Purification	70

PCR Amplification System

AmpONE [™] Taq DNA Polymerase	74
AmpONE [™] α -Taq DNA Polymerase	75
AmpONE [™] Pfu DNA Polymerase	76
AmpONE [™] HS-Taq DNA Polymerase	77
AmpONE [™] Clean Taq / Clean α -Taq DNA Polymerase	78
AmpONE [™] Taq / α -Taq / HS-Taq Premix	79
AmpMaster [™] Taq / α -Taq / HS-Taq	80
Selection Guide for DNA Polymerase	81

Index

GeneAll [®] SV Column and Filter Column Accessories	82
Ordering Information	83

GeneAll[®] Exfection[™]

For more information about products, visit www.geneall.com


1. Plasmid DNA Purification System



Hybrid-Q™ Plasmid Rapidprep	8
Exprep™ Plasmid SV mini / Midi	10
Exfection™ Plasmid LE mini / Midi	12
Exfection™ Plasmid EF Midi	12
Selection Guide for Plasmid DNA Purification	14

GeneAll[®] Hybrid-Q[™] Plasmid Rapidprep

Procedures



Component list

GeneAll[®] SV column type Qf
EzClear[™] filter
 Collection tube
 Buffer S1 - Cell Resuspension Buffer
 Buffer S2 - Cell Lysis Buffer
 Buffer G3 - Neutralization Buffer
 Buffer AW - Column Wash Buffer A
 Buffer PW - Column Wash Buffer P
 Buffer EB - Elution Buffer
 RNase A (20 mg/ml)
 Protocol Handbook

Description


GeneAll[®] Hybrid-Q[™] Plasmid Rapidprep kit provides two methods for easy and rapid preparation of plasmid DNA from the mini scale bacterial cells. Plasmid DNA can be prepared from up to 10 ml of overnight culture by conventional miniprep method with standard protocol. Alternatively, up to 3 ml of sample can be processed by rapid protocol in just 10 minutes with new patented EzClear[™] filter and simultaneous processing of multiple samples can be easily performed.

Up to 30 μ g of pure plasmid can be purified using **GeneAll[®] Hybrid-Q[™] Plasmid Rapidprep** kit and this pure plasmid DNA is ready for PCR, cloning, fluorescent sequencing, synthesis of labeled hybridization probes, cell transfection, electroporation and enzymatic restriction analysis without further manipulation.

Features and Benefits

- Spin column format
- Rapid purification with EzClear[™] filter : complete in just 10 minutes
- Stable and consistent result
- Instant use : No need of additional materials
- 30 μ g of binding capacity and high purity
- Compatible with endA⁺ strains
- No use of organic solvents
- Ready for use in fluorescent sequencing, cloning, hybridization, electroporation and other enzymatic manipulation

DNA Automated Sequencing Analysis



Plasmid DNA prepared using **GeneAll[®] Hybrid-Q[™]** delivers long and accurate (> 99 % at 700 bp) reads.
 * Sequencing analysis was performed on an ABI Prism[™], model 377, version 3.2 sequencer

For the rapid purification of high / low-copy plasmid DNA

Hybrid-Q™ Plasmid Rapidprep



Format : mini SV column type Qf,
EzClear™ filter with 2.0 ml collection tubes

Recommended sample volume (High copy) : 2 ~ 5 ml

Maximum sample volume (Low copy) : 10 ml

Maximum loading volume of EzClear™ filter : 600 µl

Maximum loading volume of spin column : 800 µl

Binding capacity : 30 µg

Recover rate : 85 ~ 95 %


Minimum elution volume : 40 µl

EzClear™ filter

New patented EzClear™ filter column facilitates the clearance of the lysate by filtration instead of tedious centrifugation which has been used widely in traditional methods.

In the rapid protocol, EzClear™ filter column is assembled with GeneAll® spin column, and this column stack makes it one-step the clearance of lysate and the binding of plasmid DNA to spin column membrane.


Stable and Consistent Result



Plasmid DNA isolated from overnight cultures of pUC18-transformed DH5αF using GeneAll® Hybrid-Q™. Each lane represents 4 µl of purified supercoiled plasmid out of 50 µl of eluates.

M : 1 Kb ladder

Enzyme Digestion



Several Kinds of plasmid DNA purified with GeneAll® Hybrid-Q™ subjected to digestion by restriction enzyme.

M : 1 Kb ladder

Lane 1 : pUC18, host DH10B

Lane 2 : pUC18, host DH10B, digested with EcoRI

Lane 3 : pUC18, host DH10B, digested with HindIII

Lane 4 : pUC18, host DH10B, digested with SmaI

Lane 5 : pQE30, host BL21

Lane 6 : pQE30, host BL21, digested with EcoRI

Lane 7 : pQE30, host BL21, digested with SalI

Lane 8 : pQE30, host BL21, digested with SmaI

Lane 9 : pACYC184, host JM109

Lane 10 : pACYC184, host JM109, digested with EcoRI

Lane 11 : pACYC184, host JM109, digested with NcoI

Lane 12 : pACYC184, host JM109, digested with PvuII


Products	Type	Size	Cat. No.
GeneAll® Hybrid-Q™ Plasmid Rapidprep	mini / spin	50	100-150
GeneAll® Hybrid-Q™ Plasmid Rapidprep	mini / spin	200	100-102

GeneAll[®] Exprep[™] Plasmid SV

Procedures

in microcentrifuges

on vacuum manifolds



Description

GeneAll[®] Exprep[™] Plasmid SV DNA Purification kit provides a rapid and convenient method for the small and medium scale preparations of plasmid DNA from bacterial cells and it is used to isolate and purify any plasmids from any E. coli strains. **GeneAll[®] Exprep[™] Plasmid SV** eliminates the need of organic solvent extraction and alcohol precipitation, allowing rapid and convenient preparation from many samples simultaneously. **GeneAll[®] Exprep[™] Plasmid SV** kit can yield up to 30 µg (mini) of highly purified plasmid DNA and it can be applicable directly for PCR, cloning, automated sequencing, synthesis of labeled hybridization probes and other enzymatic reactions without further manipulation.

Features and Benefits


- Spin or vacuum column format
- Stable and consistent result
- Instant use : No need of additional materials
- Fast and simple procedure : complete in 25 minutes
- High purity : $A_{260/280} = 1.8 \sim 2.0$
- Compatible with endA⁺ strains
- No use of organic solvents
- Ready for use in enzymatic manipulation and automated sequencing

Component list

GeneAll[®] SV column type Q
 Collection tube
 Buffer S1 - Cell Resuspension Buffer
 Buffer S2 - Cell Lysis Buffer
 Buffer S3 - Neutralization Buffer
 Buffer AW - Column Wash Buffer A
 Buffer PW - Column Wash Buffer P
 Buffer EB - Elution Buffer
 MixVu[™] Solution
 RNase A (20 mg/ml)
 Protocol Handbook

* **GeneAll[®] Midi kits** require the centrifuge which has swing bucket rotor and ability of 4,000 xg at least.

DNA Automated Sequencing Analysis



Plasmid DNA prepared using **GeneAll[®] Plasmid SV** kit delivers long, accurate (> 99 % at 700 bp) reads.

* Sequencing analysis was performed on an ABI Prism[™], model 377, version 3.2 sequencer

For the purification of plasmid DNA

Exprep™
Plasmid SV mini

Format : mini SV column type Q
with 2.0 ml collection tubes

Sample size : Up to 10 ml LB

Preparation time : 25 ~ 30 min

Typical yield : 10 ~ 30 µg

Elution volume : 30 ~ 200 µℓ

Exprep™
Plasmid SV Midi

Midi SV column type Q
with 50 ml collection tubes


Up to 100 ml LB

50 ~ 60 min

100 ~ 300 µg

400 ~ 2000 µℓ


Stable and Consistent Result



Plasmid DNA was isolated from overnight cultures of pUC18-transformed DH10B cells using **GeneAll® Exprep™** Plasmid SV kit. Each lane represents 4 µℓ of purified supercoiled plasmid DNA out of 50 µℓ of eluates.

M : 1 kb ladder

Enzyme Digestion



Several kinds of plasmid DNA purified with **GeneAll® Exprep™** Plasmid SV kit were subjected to digestion by restriction enzyme.

M : Lambda-BatPI

Lane 1 : pUC18, host INVaF

Lane 2 : pUC18, host INVaF, digested with XbaI

Lane 3 : pQE30, host JM109

Lane 4 : pQE30, host JM109, digested with SmaI

Lane 5 : pBluescript II SK (+), host XL1 blue


Lane 6 : pBluescript II SK (+), host XL1 blue, digested with XbaI

Products	Type	Size	Cat. No.
GeneAll® Exprep™ Plasmid SV	mini / spin / vacuum	50	101-150
GeneAll® Exprep™ Plasmid SV	mini / spin / vacuum	200	101-102
GeneAll® Exprep™ Plasmid SV	Midi / spin / vacuum	26	101-226
GeneAll® Exprep™ Plasmid SV	Midi / spin / vacuum	100	101-201

* Midi kit contains indicator "MixVu™" for successful alkaline lysis.

GeneAll® Exfection™ Plasmid LE / EF

Procedures



Component list

GeneAll® SV column type Qe / E
 EzClear™ Midi filter
 Collection tube
 Buffer P1
 Buffer P2
 Buffer G3 / P3
 Buffer EW1
 Buffer EW2
 Buffer ER
 Buffer EG
 Buffer EF
 MixVu™ Solution
 RNase A (20 mg/ml)
 Protocol Handbook

* GeneAll® Midi kits require the centrifuge which has a swing-bucket rotor and ability of 4,000 xg at least.

Description

GeneAll® Exfection™ plasmid LE (Low Endotoxin) and EF (Endotoxin-Free) provide simple and fast method for the purification of plasmid DNA with low endotoxin contaminants. Endotoxins (also known as lipopolysaccharides, LPS) are present in the cell membrane of gram-negative bacteria, such as Escherichia coli.


It is a common contaminant in plasmid preparations and can significantly reduce transfection efficiencies, if not removed during DNA preparations. These kits use two methods for the removal of endotoxins: advanced phase separation and endotoxin removal washing. Endotoxin levels can be reduced to 0.1 EU/μg with Exfection™ EF and to 10 EU/μg with Exfection™ LE.

Prepared plasmid DNA can be used for the transfection of most of cell lines in addition to most of molecular biological applications.

Features and Benefits

- Spin column format based on glassfiber membrane
- Convenient clearing of lysate with EzClear™ filter (Midi)
- High plasmid recoveries with high purity
- Fast preparation time and simple procedure
- High transfection efficiency in most cell-lines
- No need of additional materials
- No use of organic solvents

Transfection efficiency



pGFP-N3 prepared by the methods indicated were transfected to COS-7 (◆) and primary rat hippocampal astrocytes (■) by liposomal method. Transfection efficiencies were determined by scoring the number of green fluorescent cells 48 hours post transfection. Average transfection efficiencies are expressed as percentages relative to the efficiency obtained with DNA prepared using Exfection LE (100 %) for COS-7 and Exfection EF (100 %) for primary cells, respectively.

For the preparations of extremely pure plasmid DNA

Exfection™
Plasmid LE mini

Format : Mini column type Qe
with 2 ml collection tubes

Exfection™
Plasmid LE Midi

Midi column type E
with 50 ml collection tubes

Exfection™
Plasmid EF Midi

Midi column type E
with 50 ml collection tubes

Recommended sample volume : 5 ~ 10 ml LB

Method of lysate clearing : Centrifugation

Preparation time : < 30 min

Binding capacity : 30 µg

Endotoxin levels : < 10 EU/µg

Minimum elution volume : 50 µl

Applications : Cell transfection of most cell lines
Enzymatic modifications
Library construction
in vitro transcription/translation
High quality sequencing
Cloning
Most molecular biological experiments

50 ~ 100 ml LB

EzClear™ Midi filter

< 50 min

300 µg

< 10 EU/µg

500 µl

30 ~ 150 ml LB

EzClear™ Midi filter

< 70 min

300 µg

< 0.1 EU/µg

500 µl

In addition to LE series;
Cell transfection of
primary, sensitive and/or
suspension cell lines
Gene silencing
Microinjection

Products	Type	Size	Cat. No.
GeneAll® Exfection™ Plasmid LE	mini / spin / vacuum	50	111-150
GeneAll® Exfection™ Plasmid LE	mini / spin / vacuum	200	111-102
GeneAll® Exfection™ Plasmid LE	Midi / spin / vacuum	26	111-226
GeneAll® Exfection™ Plasmid LE	Midi / spin / vacuum	100	111-201
GeneAll® Exfection™ Plasmid EF	Midi / spin	20	121-220
GeneAll® Exfection™ Plasmid EF	Midi / spin	100	121-201


* Midi kit contains indicator "MixVu™" for successful alkaline lysis.

For Plasmid DNA Purification

Hybrid-Q™ / Exprep™ / Exfection™ Series

GeneAll® Plasmid DNA Purification Systems utilize glass microfiber membranes based on the modified alkaline lysis method. GeneAll® Hybrid-Q™ Plasmid Rapidprep with new patented EzClear™ filter provides the alternative methods for standard or rapid preparation of Plasmid DNA depending on plasmid copy number, host strain, culture medium and culture volume.

GeneAll® Exfection™ Plasmid LE (Low Endotoxin) and EF (Endotoxin-Free) provide simple and fast method for the purification of plasmid DNA with low endotoxin contaminants. Endotoxins present in the cell membrane of gram-negative bacteria are common contaminants in plasmid preparations and can significantly reduce transfection efficiencies. Exfection™ series can be used for the transfection of most cell lines through the removal of endotoxins: advanced phase separation and endotoxin removal washing.



For Plasmid DNA Purification

Hybrid-Q™ / Exprep™ / Exfection™ Series

	Hybrid-Q™ Plasmid Rapidprep*	Exprep™ Plasmid SV mini	Exfection™ Plasmid LE mini (Low Endotoxin)	Exprep™ Plasmid SV Midi	Exfection™ Plasmid LE Midi (Low Endotoxin)	Exfection™ Plasmid EF Midi (Endotoxin Free)**
	Small Scale			Large Scale		
Specifications						
Format	Spin / Vacuum	Spin / Vacuum	Spin / Vacuum	Spin / Vacuum	Spin / Vacuum	Spin
Recommended sample volume	~ 5 ml	~ 5 ml	~ 5 ml	~ 50 ml	~ 50 ml	~ 100 ml
Maximum sample volume	10 ml	10 ml	10 ml	100 ml	100 ml	150 ml
Clearing of lysate	EzClear™	Centrifuge	Centrifuge	EzClear™	EzClear™	EzClear™
Preparation time	< 10 min	< 23 min	< 30 min	< 50 min	< 50 min	< 70 min
Maximum loading volume	600 µℓ	800 µℓ	800 µℓ	15 ml	15 ml	15 ml
Binding capacity	30 µg	30 µg	30 µg	300 µg	300 µg	300 µg
The level of endotoxin	-	-	< 10 EU/µg	-	< 10 EU/µg	< 0.1 EU/µg
Recovery	85 ~ 95 %	85 ~ 95 %	80 ~ 95 %	80 ~ 95 %	85 ~ 95 %	75 ~ 90 %
Minimum elution volume	40 µℓ	50 µℓ	50 µℓ	500 µℓ	500 µℓ	500 µℓ
Applications						
Endotoxin free	-	-	-	-	-	■
Cell transfection	□	□	■	□	■	■
in vitro Transcription	■	■	■	■	■	■
Cloning	■	■	■	■	■	■
Automatic sequencing	■	■	■	■	■	■
PCR	■	■	■	■	■	■
Restriction digestion	■	■	■	■	■	■
Transformation	■	■	■	■	■	■

* Hybrid-Q™ Plasmid Rapidprep provides the alternative protocols upon plasmid copy number, host strain, culture medium, and culture volume.

** Exfection™ EF kit is suitable for the transfection of primary or sensitive cells.

• GeneAll® SV Midi/MAXI kits require the centrifuge which has a swing-bucket rotor and ability of 4,000 xg at least.

GeneAll[®] Expin[™]

For more information about products, visit www.geneall.com

2. Fragment DNA Purification System




Expin™ Gel SV	18
Expin™ PCR SV	19
Expin™ CleanUp SV	20
Expin™ Combo GP	21
Selection Guide for Fragment DNA Purification	22

GeneAll® Expin™ Gel SV

For gel extraction of
DNA fragments from agarose gel

Procedures



Component list

GeneAll® SV column type D
Collection tube
Buffer GB - Gel Extraction Buffer
Buffer NW - Column Wash Buffer N
Buffer EB - Elution Buffer
Protocol Handbook


Description

DNA extraction from agarose gel is a common technique for isolation of specific fragments from reaction mixtures. However, the conventional methods either fail to completely remove agarose, shear the DNA or result in low yields. **GeneAll® Expin™ Gel SV** kit takes advantage of glass fiber membrane to recover DNA of 80 bp to 10 kb from most grades of agarose gel in yields reaching 90 %.

Features and Benefits


- Spin or vacuum column format
- DNA extraction from standard and low-melting agarose (TAE, TBE)
- Stable and consistent result
- Rapid and convenient procedure
- Instant use : No need of additional materials
- High purity : 1.8 ~ 2.0
- Yield : 70 ~ 85 % average (80 bp ~ 10 kb)
- No use of organic solvents
- Complete removal of ethidium bromide
- pH indicator in binding buffer
- Ready for ligation, sequencing, labeling, PCR, enzyme assay and etc.

High recoveries from gels



DNA fragments resolved on 1 % agarose gel in TAE buffer.
Lane 1 ~ 5 : Before extraction with the **GeneAll® Expin™ Gel SV** kit.
Lane 6 : Pooled after extraction
* Fragment size : (up) 5.0 kb, 2.3 kb, 1.3 kb, 782 bp, 466 bp (bottom)

Extraction efficiency



Quantities of extracted 4.5 kb DNA fragment correspond to 1/5 of the DNA obtained by purification from 0.5 µg starting DNA with a recovery of 90 %. Samples were run on 1 % TAE agarose gel.
M : Lambda-BstP1
Lane 1 : Total amount before extraction (0.5 µg)
Lane 2 : 1/5 amount after extraction [90 ng compared to known amount (10 ~ 150 ng) DNA]
* Total obtained amount of DNA = 90 x 5 = 450 ng approximately (90 %)

Products	Type	Size	Cat. No.
GeneAll® Expin™ Gel SV	mini / spin / vacuum	50	102-150
GeneAll® Expin™ Gel SV	mini / spin / vacuum	200	102-102


GeneAll® Expin™ PCR SV

For the purification of DNA from PCR or other enzymatic reactions

Procedures

in microcentrifuges

on vacuum manifolds



Component list

GeneAll® SV column type D
 Collection tube
 Buffer PB - PCR Purification Buffer
 Buffer NW - Column Wash Buffer N
 Buffer EB - Elution Buffer
 Protocol Handbook

Description


GeneAll® Expin™ PCR SV kit provides fastest and easiest method for reliable purification of DNA from PCR products or other enzymatic reaction mixtures without agarose gel electrophoresis. In this kit, glass fiber membrane is used to recover DNA of 100 bp to 10 kb, which is free of primer dimers, nucleotides, enzymes and salts in yields reaching 95 %.

No organic extraction and alcohol precipitation are needed and multiple samples can be easily processed simultaneously.

Features and Benefits


- Spin or vacuum column format
- Remove PCR primers and contaminants
- Stable and consistent result
- Fast and simple : completed just in 5 minutes
- Instant use : No need of additional materials
- High purity : 1.8 ~ 2.0
- Yield : 90 ~ 95 % average (100 bp ~ 10 kb)
- No use of organic solvents
- Applied directly in ligation, automated sequencing, restriction enzyme assay, PCR, in vitro transcription, hybridization, microarray assay and other enzymatic reactions

Complete primer removal after PCR



PCR products which has several length of fragment were purified with **GeneAll® Expin™ PCR SV** kit. Enzyme, salts and small fragments such as primer dimers were effectively removed by purification. PCR product sizes : 312 bp (Lane 1, 2), 850 bp (Lane 3, 4), 1.6 kb (Lane 5, 6).
 M : phi-x174-HaeIII

Extraction efficiency




PCR products purified with **GeneAll® Expin™ PCR SV** kit was subjected to digestion with SmaI (Lane 2 ~ 5, 7 ~ 10). Lane 1, 6 represent undigested DNA.
 M : 1 kb ladder

Products	Type	Size	Cat. No.
GeneAll® Expin™ PCR SV	mini / spin / vacuum	50	103-150
GeneAll® Expin™ PCR SV	mini / spin / vacuum	200	103-102

GeneAll® Expin™ CleanUp SV

For oligonucleotide and DNA cleanup from enzymatic reactions

Procedures



Component list

- GeneAll® SV column type D
- Collection tube
- Buffer NR - DNA Clean Up Buffer
- Buffer NW - Column Wash Buffer N
- Buffer EB - Elution Buffer
- Protocol Handbook


Description

GeneAll® Expin™ CleanUp SV kit is designed for fast and simple method for purification of fragment DNA of 40 bp to 10 kb from various enzymatic reactions in just 5 minutes. Purified DNA with this kit is free of nucleotides, enzymes and salts in yields reaching 95 %, and is ready for automated sequencing, cloning, in vitro transcription, microarray and other enzymatic reactions.


Features and Benefits

- Spin or vacuum column format
- Stable and consistent result
- Fast and simple : completed just in 5 minutes
- Instant use : No need of additional materials
- High purity : 1.8 ~ 2.0
- Yield : 80 ~ 95 % (40 bp ~ 10 kb)
- No use of organic solvents
- Applied directly in ligation, automated sequencing, restriction enzyme assay, PCR, in vitro transcription, hybridization, microarray assay and other enzymatic reactions

Consistent Result from various size



DNA Automated Sequencing Analysis



Products	Type	Size	Cat. No.
GeneAll® Expin™ CleanUp SV	mini / spin / vacuum	50	113-150
GeneAll® Expin™ CleanUp SV	mini / spin / vacuum	200	113-102


GeneAll® Expin™ Combo GP

Combined kit of
Expin™ Gel SV and PCR SV

Procedures

in microcentrifuges

on vacuum manifolds



Description

GeneAll® Expin™ Combo GP kit is the combined product of Expin™ Gel SV and Expin™ PCR SV. It contains not only buffer GB required for Gel SV but also buffer PB for PCR SV, so the procedure can be chosen as user's need. No organic extraction and alcohol precipitation are needed and multiple samples can be easily processed simultaneously. Purified DNA is ready for automated sequencing, cloning, in vitro transcription, microarray and other enzymatic reaction.

Features and Benefits

- Spin and vacuum format - rapid and convenient procedure
- DNA purification from agarose gel and enzymatic reactions
- Stable and consistent result
- High yield and purity
- Instant use : No need of additional materials
- No organic extraction or alcohol precipitation
- Ready for use in cloning, automated sequencing, in vitro transcription, labeling, microarray, hybridization and other enzymatic reactions

Component list

GeneAll® SV column type D
Collection tube
Buffer GB - Gel Extraction Buffer
Buffer PB - PCR Purification Buffer
Buffer NW - Column Wash Buffer N
Buffer EB - Elution Buffer
Protocol Handbook

Recovery rates (%)

DNA size (bp)	Gel SV	PCR SV	CleanUp SV
60	39	0	63
120	71	78	80
200	76	83	84
800	84	94	94
1800	82	91	93
4300	78	85	88
8700	73	76	79


Average recovery rates of **GeneAll®** Expin™ SV kit with various sizes of DNA. 3 µg of starting sample was purified and eluted with 50 µl of buffer EB. Optional steps were not performed and SV columns were incubated for 1 minute after addition of buffer EB.

Products	Type	Size	Cat. No.
GeneAll® Expin™ Combo GP	mini / spin / vacuum	50	112-150
GeneAll® Expin™ Combo GP	mini / spin / vacuum	200	112-102

For Fragment DNA Purification

Expin™ Series

GeneAll® Expin™ series provide reliable and fast methods for the purification of fragment DNA from agarose gel and PCR or enzymatic reaction mixtures. Expin™ Gel SV takes advantage of glass fiber membrane to recover DNA of 80 bp to 10 kb from most grades of agarose gel in yields reaching 90 %. Expin™ PCR SV is used to recover DNA of 100 bp to 10 kb from PCR or enzymatic reaction mixtures and very effective to the removal of PCR primer dimer. Expin™ Combo GP kit is the combined product of Expin™ Gel SV and Expin™ PCR SV. Expin™ CleanUp SV is designed for fast and simple method for purification of fragment DNA of 40 bp to 10 kb from various enzymatic reactions in just 5 minutes.



* GeneAll® Expin™ SV series consist of Gel, PCR and CleanUp SV kit. Each kit is optimized for efficient recovery of DNA and removal of contaminants in each specific application.

	Expin™ PCR SV	Expin™ Gel SV	Expin™ CleanUp SV	Expin™ Combo GP*
Specifications				
Format	Spin / Vacuum	Spin / Vacuum	Spin / Vacuum	Spin / Vacuum
Starting material	100 µl PCR reactions	200 mg Gel slice	50 µl Enzyme reactions	100 µl PCR reactions 200 mg Gel slice
Fragment DNA size	100 bp ~ 10 kb	80 bp ~ 10 kb	40 bp ~ 10 kb	80 bp ~ 10 kb
Recovery	90 ~ 95 %	70 ~ 85 %	80 ~ 95 %	70 ~ 95 %
Maximum binding capacity	10 µg	10 µg	10 µg	10 µg
Preparation time	< 5 min	< 15 min	< 5 min	< 5 min ~ 15 min
Applications				
PCR cleanup	■	-	□	■
Gel extraction	-	■	-	■
Nucleotide removal	□	□	■	□

■ Recommended / □ Suitable but not optimized

* Expin™ Combo GP kit is the combined product of Expin™ Gel SV and Expin™ PCR SV.

3. Genomic DNA Purification System




Exgene™ Tissue SV (<i>plus!</i>)	24
Exgene™ Blood SV	26
Exgene™ Clinic SV	28
Exgene™ Cell SV	30
Exgene™ Plant SV	32
Exgene™ Soil DNA mini	34
Exgene™ Genomic DNA micro	36
Exgene™ GMO SV	38
GenEx™ Blood / Cell / Tissue	40
GenEx™ Plant (<i>plus!</i>)	42
DirEx™	44
Selection Guide for Genomic DNA Purification	46

GeneAll® Exgene™ Tissue SV (plus!)

Procedures

in microcentrifuges

on vacuum manifolds



Component list

GeneAll® SV column type G
 Collection tube
 Additional collection tube
 Buffer RL - RBC Lysis Buffer (plus! only)
 Buffer TL - Tissue Lysis Buffer
 Buffer TB - Tissue Binding Buffer
 Buffer BW - Column Wash Buffer B
 Buffer TW - Column Wash Buffer T
 Buffer AE - Elution Buffer
 Proteinase K (20 mg/ml)
 Protocol Handbook

* GeneAll® Midi / MAXI kits require the centrifuge which has swing bucket rotor and ability of 4,000 xg at least.


Description

GeneAll® Exgene™ Tissue SV kit provides a simple and rapid method for the isolation of total DNA from animal tissues and cultured cells. This kit can process 25 mg (mini) of wet tissue and yields up to 50 µg (mini) depending on the type of sample used. Specially formulated buffer system minimize RNA copurified with DNA without RNase A treatment. If RNA-free genomic DNA is required, RNase A can be treated with. No organic extraction and alcohol precipitation are needed and multiple samples can be easily processed simultaneously. Exgene™ Tissue SV plus! offers additional material and method for DNA purification from whole blood.

Features and Benefits

- Spin or vacuum column format
- Accurate and consistent DNA extraction from animal tissues, cultured cell line and whole blood (plus! only)
- Instant use : No need of additional materials
- Simple and safe procedure
- High purity : 1.8 ~ 2.0
- No use of organic solvents
- Ready for use in PCR, southern blotting, AFLP, RFLP, RAPD and other enzymatic reactions

DNA Extraction from various samples



DNA purification using GeneAll® Exgene™ Tissue SV kit. DNA from several kinds of animal tissue was prepared. Elution was performed with 100 µl of Buffer AE. 8 µl of eluates was resolved on 0.8 % agarose gel.

For the isolation of gDNA from tissues, cells and whole blood (plus!)

**Exgene™
Tissue SV (plus!)
mini**



Format : mini SV column type G
with 2.0 ml collection tubes

Sample size : Up to 25 mg tissue

Preparation time : 25 ~ 220 min

Typical yield : 5 ~ 50 µg

Elution volume : 30 ~ 400 µℓ

**Exgene™
Tissue SV (plus!)
Midi**



Midi SV column type G
with 15 ml collection tubes

Up to 100 mg tissue

40 ~ 250 min

20 ~ 150 µg

200 ~ 600 µℓ

**Exgene™
Tissue SV (plus!)
MAXI**



MAXI SV column type G
with 50 ml collection tubes


Up to 250 mg tissue

40 ~ 250 min

80 ~ 400 µg

400 ~ 2000 µℓ

Restriction Enzyme Assay




Genomic DNA purified from various rat samples using **GeneAll® Exgene™ Tissue SV kit** was partially digested with **EcoRI** (Lane 2 ~ 3, 5 ~ 6, 8 ~ 9).

Lane 1, 4, 7 represent undigested DNA.

M : 1 kb ladder

PCR Amplification



PCR reaction was performed with purified DNA using **GeneAll® Exgene™ Tissue SV kit**. Template was isolated from rat liver (Lane 1 ~ 2, 7 ~ 8), spleen (Lane 3 ~ 4, 9 ~ 10) and kidney (Lane 5 ~ 6, 11 ~ 12).

M : 1 kb ladder


Products	Type	Size	Cat. No.
GeneAll® Exgene™ Tissue SV (plus!)	mini / spin / vacuum	100	104(9)-101
GeneAll® Exgene™ Tissue SV (plus!)	mini / spin / vacuum	250	104(9)-152
GeneAll® Exgene™ Tissue SV (plus!)	Midi / spin / vacuum	26	104(9)-226
GeneAll® Exgene™ Tissue SV (plus!)	Midi / spin / vacuum	100	104(9)-201
GeneAll® Exgene™ Tissue SV (plus!)	MAXI / spin / vacuum	10	104(9)-310
GeneAll® Exgene™ Tissue SV (plus!)	MAXI / spin / vacuum	26	104(9)-326

GeneAll® Exgene™ Blood SV

Procedures

in microcentrifuges

on vacuum manifolds



Component list

GeneAll® SV column type G
Collection tube
Additional collection tube
Buffer BL - Blood Lysis Buffer
Buffer BW - Column Wash Buffer B
Buffer TW - Column Wash Buffer T
Buffer AE - Elution Buffer
Proteinase K (20 mg/ml)
Protocol Handbook

* GeneAll® Midi / MAXI kits require the centrifuge which has swing bucket rotor and ability of 4,000 xg at least.


Description

GeneAll® Exgene™ Blood SV kit provides a simple and rapid method for the isolation of total DNA from fresh or frozen whole blood, buffy coat, serum, plasma, virus and cultured cells. Purification procedure is so simple and optimized to simultaneous processing of multiple samples. Exgene™ Blood SV yields pure DNA ready for direct PCR in just 20 minutes (mini) and 1 hour (Midi / MAXI). There is no need phenol extraction or alcohol precipitation.

Features and Benefits

- Spin or vacuum column format
- Accurate and consistent DNA extraction from whole blood, buffy coat, serum, plasma, cultured cells
- Instant use : No need of additional materials
- Fast, safe and simple procedure completed in 20 minutes (mini), 1 hour (Midi, MAXI)
- High purity : 1.8 ~ 2.0
- No use of organic solvents
- Ready for use in PCR, southern blotting and other enzymatic reactions

DNA Extraction from various samples



Total DNA was isolated from 200 μl of whole blood of various species using GeneAll® Blood SV mini kit. Each lane represents 8 μl of 100 μl eluates.

M : Lambda-HindIII

For the purification of gDNA from blood and its derivatives

Exgene™
Blood SV mini

Format : mini SV column type G
with 2.0 ml collection tubes

Sample size : Up to 200 μl
(Whole blood)

Preparation time : 20 ~ 30 min

Typical yield : 4 ~ 20 μg

Elution volume : 30 ~ 400 μl

Exgene™
Blood SV Midi

Midi SV column type G
with 15 ml collection tubes

1 ~ 2 ml

40 ~ 55 min

20 ~ 80 μg

200 ~ 600 μl

Exgene™
Blood SV MAXI

MAXI SV column type G
with 50 ml collection tubes


3 ~ 10 ml

40 ~ 55 min

80 ~ 400 μg

400 ~ 2000 μl


Restriction Enzyme Assay



Genomic DNA purified from various rat blood samples using **GeneAll® Exgene™ Blood SV mini** kit was partially digested with **EcoRI** (Lane 2 ~ 3, 5 ~ 6, 8 ~ 9). Lane 1, 4, 7 represent undigested DNA.

M : 1 kb ladder

PCR Amplification



PCR reaction was performed with purified DNA using **GeneAll® Exgene™ Blood SV** kit as template. Each lane 1, 2 and 3 corresponds to the template of each PCR product (Lane 4, 5, 6). Template DNA was isolated from whole blood of rat (SD) and the exon region of GAPDH gene was amplified with Taq polymerase.


Products	Type	Size	Cat. No.
GeneAll® Exgene™ Blood SV	mini / spin / vacuum	100	105-101
GeneAll® Exgene™ Blood SV	mini / spin / vacuum	250	105-152
GeneAll® Exgene™ Blood SV	Midi / spin / vacuum	26	105-226
GeneAll® Exgene™ Blood SV	Midi / spin / vacuum	100	105-201
GeneAll® Exgene™ Blood SV	MAXI / spin / vacuum	10	105-310
GeneAll® Exgene™ Blood SV	MAXI / spin / vacuum	26	105-326

GeneAll[®] Exgene[™] Clinic SV

Procedures

in microcentrifuges

on vacuum manifolds



Component list

GeneAll[®] SV column type G

Collection tube

Buffer CL - Cell Lysis Buffer

Buffer BL - Blood Lysis Buffer

Buffer BW - Column Wash Buffer B

Buffer TW - Column Wash Buffer T

Buffer AE - Elution Buffer

Proteinase K (20 mg/ml)

Protocol Handbook

* **GeneAll[®] Midi / MAXI kits** require the centrifuge which has swing bucket rotor and ability of 4,000 xg at least.


Description

GeneAll[®] Exgene[™] Clinic SV kit provides an easy and fast method for the isolation of total DNA such as genomic, mitochondrial, bacterial, parasite or viral DNA from various clinical sample including tissues, whole blood and body fluids. The purified DNA is suitable for PCR, blotting, RFLP, RAPD, AFLP and etc.

Features and Benefits

- Easy and fast purification of high-quality DNA
- Spin and vacuum format
- Instant use : No need of additional materials
- No organic extraction or alcohol precipitation
- Consistent and high yields
- High purity : 1.8 ~ 2.0
- Ready for use in PCR, Southern blotting, genotyping and etc.

Consistent Result from various sample



Total DNA purified from various sample tissue using **GeneAll[®] Exgene[™] Clinic SV mini** is resolved on 0.8 % agarose gel.
M : Lambda-HindIII

For the isolation of gDNA from clinical tissues including whole blood

Exgene™
Clinic SV mini

Format : mini SV column type G
with 2.0 ml collection tubes

Sample size : Up to 20 mg*

Preparation time : 25 ~ 220 min*

Typical yield : 5 ~ 50 µg

Elution volume : 30 ~ 400 µℓ

Exgene™
Clinic SV Midi

Midi SV column type G
with 15 ml collection tubes

Up to 100 mg*

40 ~ 250 min*

20 ~ 150 µg

200 ~ 600 µℓ

Exgene™
Clinic SV MAXI

MAXI SV column type G
with 50 ml collection tubes

Up to 250 mg*

40 ~ 250 min*

80 ~ 400 µg

400 ~ 2000 µℓ

* Depending on the sample used.


Products	Type	Size	Cat. No.
GeneAll® Expin™ Clinic SV	mini / spin / vacuum	100	108-101
GeneAll® Expin™ Clinic SV	mini / spin / vacuum	250	108-152
GeneAll® Expin™ Clinic SV	Midi / spin / vacuum	26	108-226
GeneAll® Expin™ Clinic SV	Midi / spin / vacuum	100	108-201
GeneAll® Expin™ Clinic SV	MAXI / spin / vacuum	10	108-310
GeneAll® Expin™ Clinic SV	MAXI / spin / vacuum	26	108-326

GeneAll[®] Exgene[™] Cell SV

Procedures

in microcentrifuges

on vacuum manifolds



Component list

GeneAll[®] SV column type G
 Collection tube
 Buffer GP - Enzyme Reaction Buffer
 Buffer YL - Lyticase Reaction Buffer
 Buffer CL - Cell Lysis Buffer
 Buffer BL - Blood Lysis Buffer
 Buffer BW - Column Wash Buffer B
 Buffer TW - Column Wash Buffer T
 Buffer AE - Elution Buffer
 Proteinase K (20 mg/ml)
 Protocol Handbook

* **GeneAll[®] MAXI Kits** require the centrifuge which has swing bucket rotor and ability of 4,000 xg at least.


Description

GeneAll[®] Cell SV kit provides a rapid and simple method for the purification of total DNA from a wide range of organism including bacterial cells, yeast, cultured cells, whole blood and blood derivatives. Up to 2×10^9 bacterial cells, 5×10^6 cultured cells or 3 ml of yeast cultures may yield 10 ~ 40 μ g of DNA typically. The pure DNA can be acquired in just 30 minutes and this can be directly used in various applications such as PCR, Southern blotting and other enzymatic reactions.

Features and Benefits

- Spin or vacuum column format
- Accurate and consistent DNA extraction from gram positive or negative bacteria, cultured cell, yeast and various biological samples
- High purity : 1.8 ~ 2.0
- Simple and safe procedure
- No use of organic solvents
- Ready for use in PCR, Southern blotting, AFLP, RFLP, RAPD and other enzymatic reactions

Consistent Result from various samples



Genomic DNA prepared from a several species of cells using **GeneAll[®] Exgene[™] Cell SV kit**. 5 μ l out of 100 μ l eluate was resolved on 0.8 % agarose gel.

M : Lambda-HindIII

For the isolation of gDNA from cultured cell, yeast
gram positive/negative bacteria and etc.

Exgene™ Cell SV mini



Format : mini SV column type G
with 2.0 ml collection tubes

Sample size : Up to 2×10^9 bacterial cells
Up to 5×10^7 yeast cells

Preparation time : 30 ~ 120 min*

Typical yield : 10 ~ 50 μg

Elution volume : 30 ~ 400 μl

Exgene™ Cell SV MAXI



MAXI SV column type G
with 50 ml collection tubes

Up to 2×10^{10} bacterial cells
Up to 5×10^8 yeast cells


60 ~ 240 min*

80 ~ 400 μg

400 ~ 2000 μl

* Depending on the sample used.


Restriction Enzyme Assay



Genomic DNA purified from E. coli DH5 α F and JB6 samples using GeneAll® Exgene™ Cell SV kit was partially digested with BamHI (Lane 2 ~ 4, 6 ~ 8). Lane 1, 5 represent undigested DNA.

M : Lambda-HindIII

PCR Amplification



PCR reaction was performed with purified DNA using GeneAll® Exgene™ Cell SV kit. Template DNA was isolated from E. coli DH5 α F (Lane 1, 2, 3).

PCR reaction was performed with genomic DNA purified from DH5 α F using GeneAll® Cell SV kit.

M : Lambda-HindIII


Products	Type	Size	Cat. No.
GeneAll® Exgene™ Cell SV	mini / spin / vacuum	100	106-101
GeneAll® Exgene™ Cell SV	mini / spin / vacuum	250	106-152
GeneAll® Exgene™ Cell SV	MAXI / spin / vacuum	10	106-310
GeneAll® Exgene™ Cell SV	MAXI / spin / vacuum	26	106-326

GeneAll[®] Exgene[™] Plant SV

Procedures

in microcentrifuges

on vacuum manifolds



Component list

GeneAll[®] SV column type G

Collection tube
 EzSep[™] filter column
 Collection tubes
 Buffer PL - Plant Lysis Buffer
 Buffer PD - Precipitation Buffer
 Buffer BD - Binding Buffer
 Buffer CW - Column Wash Buffer C
 Buffer AE - Elution Buffer
 RNaseA (100 mg/ml)
 Protocol Handbook

* **GeneAll[®] Midi / MAXI Kits** require the centrifuge which has swing bucket rotor and ability of 4,000 xg at least.


Description

GeneAll[®] Exgene[™] Plant SV kit provides a simple and easy method for the small, medium and large scale purification of total DNA from various plant tissues. With EzSep[™] filter and **GeneAll[®]** column type G, the procedure can be done in just 40 minutes (mini), yielding a pure genomic DNA suitable for various downstream applications without further manipulation. Up to 100 mg, 400 mg of plant tissue can be processed with Exgene[™] Plant SV mini, Midi and MAXI, respectively. Exgene[™] Plant SV procedure eliminates the need of organic solvent extraction and alcohol precipitation, allowing safe and fast preparation of many samples simultaneously. Purified total DNA can be directly applicable in conventional PCR, real-time PCR, Southern blotting, SNP genotyping, RFLP, AFLP and RAPD.

Features and Benefits

- Spin or vacuum column format
- Stable and consistent DNA extraction from plant cells, tissues and fungi
- Instant use : No need of additional materials
- Perfect removal of second metabolites such as polyphenols and polysaccharides
- Simple procedure by the use of EzSep[™] filter column
- No use of organic solvents
- Ready for use in PCR, Southern blotting, AFLP, RFLP, RAPD and other enzymatic reactions

Comparison of DNA Extraction with supplier Q



Genomic DNA was extracted from each 100 mg of various samples and analyzed on 0.8 % agarose gel. To compare with supplier Q, same kind and amount of each plant samples were subjected to extraction.
 Q : supplier Q, G : **GeneAll[®] Exgene[™] Plant SV** kit.
 M : 1 kb ladder

For the isolation of gDNA from plant cells and tissues

Exgene™
Plant SV miniExgene™
Plant SV MidiExgene™
Plant SV MAXI

Format : mini SV column type G
with 2.0 ml collection tubes

Midi SV column type G
with 15 ml collection tubes

MAXI SV column type G
with 50 ml collection tubes

Sample size : Up to 100 mg wet (30 mg dry)

Up to 400 mg wet (100 mg dry)

Up to 1 g wet (300 mg dry)

Preparation time : < 40 min

< 1 hour

< 1 hour

Typical yield : 4 ~ 40 µg

10 ~ 150 µg


40 ~ 300 µg

Elution volume : 30 ~ 400 µl

200 ~ 600 µl

400 ~ 2000 µl


Restriction Enzyme Assay



Genomic DNA purified from various plant samples by GeneAll® Exgene™ Plant SV kit was subjected to partial digestion with HindIII (Lane 2, 4, 6, 8). Lane 1, 3, 5, 7 represent undigested DNA.

M : 1 kb ladder

PCR Amplification




PCR reaction was performed with purified DNA using GeneAll® Exgene™ Plant SV kit. Two primer sets were used : trnL region (left lanes) and large subunit rDNA region on plasmid

Products	Type	Size	Cat. No.
GeneAll® Exgene™ Plant SV	mini / spin / vacuum	100	117-101
GeneAll® Exgene™ Plant SV	mini / spin / vacuum	250	117-152
GeneAll® Exgene™ Plant SV	Midi / spin / vacuum	26	117-226
GeneAll® Exgene™ Plant SV	Midi / spin / vacuum	100	117-201
GeneAll® Exgene™ Plant SV	MAXI / spin / vacuum	10	117-310
GeneAll® Exgene™ Plant SV	MAXI / spin / vacuum	26	117-326

GeneAll® Exgene™ Soil DNA mini

Procedures



Component list

GeneAll® Column type G
(with collection tube)
1.5 ml collection tube
2.0 ml collection tube
Buffer SL
Buffer RH
Buffer PD
Buffer TB
Buffer NW
Buffer EB
Powerbead™ tube
Protocol Handbook

Description

GeneAll® Exgene™ Soil DNA mini provides a convenient method for the isolation of total DNA from soil samples. This kit utilizes the powerful beads, the optimized buffer system and the advanced silica binding technology to purify nucleic acid suitable for many applications. These complex systems of this kit can deal with a number of different types of samples in the soil including plant tissues, bacteria, fungi spores and others. Also, it removes a humic acid contents and other PCR inhibitors from various soil samples efficiently. The humic acid contents which is a sort of brownish colour, is a critical factor for soil treating experiments. If remained in eluate, this can have a negative effect on the DNA downstream applications.

GeneAll® Exgene™ Soil DNA mini provides a tube including powerful beads for strong pulverization. Soil samples are placed in this tube with lysis buffer, buffer SL, and crushed by bead-beater or vortex. After centrifugation, supernatant is mixed with precipitation buffer, buffer RH and buffer PD, to precipitate humic acid and protein. Then, the separated DNA part, supernatant, blend into the binding buffer, buffer TB, and DNA is bound on the silica membrane through centrifugation. Following washing step with buffer NW, the bound DNA is eluted by buffer EB. Purified DNA can be directly applicable in conventional PCR, restriction analysis, electrophoresis and any other downstream applications.

Features and Benefits

- Glassfiber membrane technology
- Sample size : Up to 500 mg
- Easy and fast purification of high-quality DNA
- Preparation time : ~ 25 minutes
- Efficient lysis step using Powerbead™ tube
- Perfect removal of humic acid
- Stable and consistent yield
- No organic extraction or alcohol precipitation
- High purity : ready for the conventional and real-time PCR

For the isolation of gDNA from soil samples

Exgene™
Soil DNA mini



Format : mini SV column type G with
2.0 ml collection tubes

Sample size : Up to 500 mg

Preparation time : ~ 25 min

Elution volume : 30 ~ 200 μl

Comparative Genomic DNA purification result



gDNA isolated from various soil samples with Exgene™ Soil DNA mini (Fig. I) vs
supplier A (Fig. II) (used vortex homogenization method)

Lane M : Lambda-HindIII


Lane 1 : Soil under cherry blossom

Lane 2 : Soil of onion patch

Lane 3 : Soil of cabbage patch

Lane 4 : Mud

PCR result



gDNA was purified from various soil samples using Exgene™ Soil DNA mini.
And then, the 16s rRNA was amplified by PCR and confirmed by electrophoresis.

Lane M : 100 bp ladder

Lane 1 : Pot soil

Lane 2 : Soil under cherry blossom A

Lane 3 : Soil of cabbage patch A

Lane 4 : Soil under cherry blossom B

Lane 5 : Soil of cabbage patch B

Lane 6 : Soil under cherry blossom C


Lane 7 : Soil of cabbage patch C

Lane 8 : Negative

Products	Type	Size	Cat. No.
GeneAll® Exgene™ Soil DNA mini	mini / spin	50	114-150

GeneAll[®] Exgene[™] Genomic DNA micro

Procedures



Component list

GeneAll[®] Micro column type S
 Collection tube
 Buffer CL - Lysis Buffer I
 Buffer BL - Lysis Buffer II
 Buffer BW - Column Wash Buffer B
 Buffer TW - Column Wash Buffer T
 Buffer AE - Elution Buffer
 Carrier RNA
 Proteinase K
 PK Storage Buffer
 Protocol Handbook

Description


GeneAll[®] Exgene[™] Genomic DNA micro kit provides fast and easy methods for the micro scale purification of total (genomic and mitochondrial) DNA from various biological samples. Purified DNA can be used directly for PCR, quantitative PCR, genotyping such as STR analysis and other downstream applications.

Exgene[™] Genomic DNA micro utilizes the advanced silica-binding technology to purify total DNA sufficiently pure for many applications. Various samples are lysed in optimized buffer containing detergents and lytic enzyme. Under high salt condition, DNA in the lysate bind to silica membrane and impurities pass through membrane into a collection tube. The membranes are washed with a series of alcohol-containing buffer to remove any traces of proteins, cellular debris and salts. Finally pure DNA is released into a clean collection tube with deionized water or low ionic strength buffer.

Features and Benefits

- Spin column format
- Apply to trace of sample : use of micro column
- Simple and safe procedure
- Stable and consistent result
- Instant use : No need of additional materials
- No use of organic solvents
- High yield and purity
- Various protocol for forensic sample : stain, chewing gum, cigarette butts, tooth brush
- Ready for use in general PCR, qPCR, genotyping such as STR analysis and other downstream applications

Stable and Consistent Result



Automatic sequencing data of 1 kb PCR products of extracted genomic DNA by **GeneAll[®] Exgene[™] Genomic DNA micro** kit. Sequencing was performed on an ABI3730XL (96-capillary) DNA sequencer using an internal primers.

For the isolation of total DNA from micro-scale biological samples

Exgene™ Genomic DNA micro




Format : Micro column type S
with 2.0 ml collection tubes

Sample size : Up to 100 μl whole blood


Preparation time : After lysis, 20 min

Elution volume : 20 ~ 50 μl

Real-Time PCR Amplification




Real-time PCR was performed with purified DNA **GeneAll® Exgene™** Genomic DNA micro kit. The DNA was extracted from whole blood, stains, swab and hair root (Panel 1), nail clippings, chewing gum, tooth brush and urine (Panel 2). Real-time PCR was carried out with human GAPDH primer sets and detected by SYBR® Green reagent.



DNA extraction from various biological samples using **GeneAll® Exgene™** Genomic DNA micro kit or a kit from Supplier A. Real-time PCR was carried out human GAPDH primer sets or mitochondria hypervariable region I primer sets and detected by SYBR® Green reagent.

PCR amplification



PCR reaction was performed with purified DNA using **GeneAll® Exgene™** Genomic DNA micro kit. Template was isolated from whole blood (Lane 1), dried blood spot (Lane 2), hair root (Lane 3), chewing gum (Lane 4), animal tissue (Lane 5), urine (Lane 6), bone (Lane 7) and hair shaft (Lane 8). M: 1 kb ladder


Products	Type	Size	Cat. No.
GeneAll® Exgene™ Genomic DNA micro	mini / spin	50	118-050

GeneAll[®] Exgene[™] GMO SV

Procedures

in microcentrifuges

on vacuum manifolds



Component list

GeneAll[®] SV column type S
 Collection tube
 Buffer GL - GMO Lysis Buffer
 Buffer MP - GMO Protein Precipitation Buffer
 Buffer MB - GMO Binding Buffer
 Buffer GW - Column Wash Buffer G
 Buffer EB - Elution Buffer
 Protocol Handbook


Description

GeneAll[®] GMO SV kit provides fast and convenient methods for the isolation of genomic DNA from the crops or its processed food such as soybean, maize, rice, barely, soybean paste, bean curd, snacks and etc. Up to 300 mg of crop can be processed with this kit in just an hour, depending on the lysis conditions. Phenol / chloroform extraction, alcohol precipitation and cold conditions are not required. Without any further manipulations, purified DNA is ready for the conventional and real-time PCR for the detection of GMO (genetically modified organism).

Features and Benefits

- Spin or vacuum column format
- Accurate and consistent DNA extraction from 300 mg of various grains or its processed food
- Fast and simple procedure
- No organic extraction or alcohol precipitation
- High purity : ready for the conventional and real-time PCR

Consistent Result from various sample



Total DNA was prepared from 150 mg of soybean (Lane 1 ~ 3) and maize (Lane 4 ~ 6) using **GeneAll[®] Exgene[™] GMO SV kit**. Each 8 μ l out of 50 μ l of elute was resolved on 0.8 % agarose gel.
 M : Lambda-HindIII

For the isolation of gDNA from various crops and its processed foods

Exgene™ GMO SV mini



Format : mini SV column type S
with 2.0 ml collection tubes


Sample size : Up to 300 mg grains

Preparation time : 40 ~ 50 min

Typical yield : 6 ~ 30 µg

Elution volume : 30 ~ 200 µℓ

GMO gene Screening test




GMO gene of soybean was detected using previous DNA as template. PCR reaction was carried out doubly (Lane 1 ~ 2, 4 ~ 5, 7 ~ 8, respectively) and its template DNA corresponds to previous DNA (Lane 1 ~ 3). Lane 3, 6 and 9 represent negative controls.

Products	Type	Size	Cat. No.
GeneAll® Exgene™ GMO SV	mini / spin / vacuum	50	107-150
GeneAll® Exgene™ GMO SV	mini / spin / vacuum	200	107-102

GeneAll[®] GenEx[™] Blood / Cell / Tissue

Procedures



Description

GenEx[™] Series provide convenient methods for the isolation of total DNA from various biological samples without use of toxic chemical such as phenol or chloroform. These kits utilize the specially formulated buffer system in order to process the sample scalably and obtain the almost intact size of genomic DNA. Extracted genomic DNA can be applied directly to PCR, southern blotting and restriction enzyme assay and other downstream applications.


GenEx[™] Series can be used for:

GenEx[™] Blood - Whole blood or blood derivatives
 GenEx[™] Cell - Cultured cells or gram negative bacteria
 GenEx[™] Tissue - Animal tissues

Features and Benefits

- Specially formulated buffer system
- DNA preparation from diverse sample; whole blood, cultured cell, yeast, bacteria, animal tissue and etc.
- Recovery of very high molecular weight DNA
- Rescalable preparation depending on sample amount
- No organic extraction
- High purity : ready for PCR, southern blotting and other downstream applications

DNA Extraction from various samples



Genomic DNA prepared from several kinds of organism using GenEx[™] Genomic DNA Isolation kit. 5 μ l of eluate from each sample was resolved on 0.7% agarose gel.

Component list

Buffer RL - RBC Lysis Buffer
 Buffer AL - Cell Lysis Buffer
 Buffer PP - Protein Precipitation Buffer
 Buffer RE - DNA Rehydration Buffer
 RNase A (20 mg/ml)
 Proteinase K
 PK storage Buffer
 Protocol Handbook


For the isolation of gDNA from whole blood, cultured cell, animal tissue and etc.

DNA Yields from various starting materials

Materials	Species	Amount	Yields of DNA
Whole blood*	Human	300 μl	5 ~ 15 μg
		3 ml	80 ~ 150 μg
		10 ml	250 ~ 500 μg
	Mouse	300 μl	6 ~ 7 μg
Buffy coat*	Human	150 ~ 250 μl	50 ~ 150 μg
Body fluids	Human	50 μl	0.1 ~ 2.5 μg
Cultured cell lines	CHO	2 x 10 ⁶ cells	14 ~ 16 μg
	RAW264.7	2 x 10 ⁶ cells	16 ~ 17 μg
	COS	1.5 x 10 ⁶ cells	9 ~ 12 μg
	K562	3 x 10 ⁶ cells	15 ~ 30 μg
	NIH3T3	2 x 10 ⁶ cells	9 ~ 13 μg
	PC12	8 x 10 ⁶ cells	5 ~ 8 μg
Animal tissue	Mouse Liver	10 mg	20 ~ 25 μg
	Mouse Pancreas	10 mg	70 ~ 75 μg
	Mouse Heart	10 mg	2 ~ 4 μg
	Mouse Tail	1 cm of tail tip	15 ~ 30 μg
Gram (-) bacteria	E.coli / JM109	2 x 10 ⁹ cells	18 ~ 25 μg
	E.cloacae	6 x 10 ⁹ cells	20 ~ 26 μg

* Yield depends on the quantity of white blood cells present.


Restriction Enzyme Assay



Genomic DNA purified from various organism samples using GenEx™ Genomic Isolation kit was partially digested with HindIII (Lane 2, 4, 6, 8, 10, and 12). Lane 1, 3, 5, 7, 9 and 11 represent undigested DNA.

M : Lambda-HindIII

PCR Amplification



PCR reaction was performed with purified DNA using GenEx™ Genomic Isolation kit. Template DNA was isolated from Tobacco (Lane 1), BL21 (Lane 2), DH5α (Lane 3), Liver (Lane 4), Kidney (Lane 5), Whole blood (Lane 6) and JB6 (Lane 7).

M : 1 kb ladder


Products	Type	Size	Cat. No.
GeneAll® GenEx™ Blood	Sx / Solution	Sx / 100	220-101
GeneAll® GenEx™ Blood	Sx / Solution	Sx / 500	220-105
GeneAll® GenEx™ Blood	Lx / Solution	Lx / 100	220-301
GeneAll® GenEx™ Cell	Sx / Solution	Sx / 100	221-101
GeneAll® GenEx™ Cell	Sx / Solution	Sx / 500	221-105
GeneAll® GenEx™ Cell	Lx / Solution	Lx / 100	221-301
GeneAll® GenEx™ Tissue	Sx / Solution	Sx / 100	222-101
GeneAll® GenEx™ Tissue	Sx / Solution	Sx / 500	222-105
GeneAll® GenEx™ Tissue	Lx / Solution	Lx / 100	222-301

Sx On the basis of DNA purification from 300 μl whole blood, 2 x 10⁶ cells or 10 mg animal tissue

Lx On the basis of DNA purification from 10 ml whole blood, 1 x 10⁹ cells or 100 mg animal tissue

GeneAll[®] GenEx[™] Plant (plus!)

Procedures



Description


GeneAll[®] GenEx[™] Plant kit provides an easy and convenient method for the isolation of total DNA from various plant samples without use of toxic chemical such as phenol or chloroform. This kit has a specially formulated solution format and enables the scalable preparation of almost intact size DNA. Especially when purifying DNA from plant, the removal of secondary metabolites is very important because contamination of these impurities can lead to inhibition of downstream application. The optimized buffer system adopted in this kit can facilitate the removal of contaminants, such as second metabolites and other impurities. Purified DNA can be applied directly to PCR, blotting, restriction enzyme assay and other downstream applications.

GeneAll[®] GenEx[™] Plant plus! kit has an additional feature, EzSep[™] filter column. With certain plant samples, it is very difficult to separate cleared supernatant from pelletal debris at a protein precipitation stage. This problem also appears often when large starting sample and it may be due to low density of debris and/or low centrifugal force with conventional centrifuge. EzSep[™] filter column included in the plus! kit is the device to solve this problem and moreover it decreases the preparation time also.

Features and Benefits

- Specially formulated buffer system
- DNA preparation from various plant samples
- Recovery of very high molecular weight DNA
- Rescalable preparation depending on sample amount
- No organic extraction
- High purity : ready for PCR, Southern blotting and other downstream applications
- Simple separation of supernatant by EzSep[™] filter (plus! only)

Result from various samples



Total DNA prepared from various plant leaves using GenEx[™] Plant kit. Each sample is extracted from 100 mg of tissue approximately. And 4 uls of purified DNA were resolved on 1.0 % agarose gel.

M : 1 kb DNA ladder


Component list

Buffer PL
Buffer PP
Buffer RE
RNase A (100 mg/ml)
EzSep[™] filter (plus! only)
Protocol Handbook

For the isolation of total DNA from various plant samples

PCR Amplification

M 1 2 3 4 5 6 7 8 9 10 11 12 13 14 N1 N2 M




PCR was performed with total DNA purified from various samples using GenEx™ Plant as template. The primer set is for a 297 bp fragment of a highly conserved region of chloroplast DNA. PCR products were resolved on 1.2% agarose gel.

Lane M : 100 bp DNA ladder

Lane 1 : Zelkova Lane 5 : Taxodium Lane 9 : Prunus Lane 13 : Citrus
Lane 2 : Lingustrum Lane 6 : Nicotiana Lane 10 : Rosa Lane 14 : Actinidia
Lane 3 : Arabidopsis Lane 7 : Gingko Lane 11 : Solanum Lane N1 : Negative control 1 – no template.
Lane 4 : Diospyros Lane 8 : Lactuca Lane 12 : Capsicum Lane N2 : Negative control 2 – E. coli gDNA as template.

Restriction Enzyme Assay

M 1 2 3 4 5 6



Total DNA(Lane 1, 3, 5) purified from the leaves of several species using GenEx™ Plant was subjected to restricted digestion(Lane 2, 4, 6) by HindIII.

Lane M : 1 kb DNA ladder

Lane 1 : Zelkova
Lane 3 : Taxodium
Lane 5 : Nicotiana

Products	Type	Size	Cat. No.
GeneAll® GenEx™ Plant	Solution	Sx	227-101
GeneAll® GenEx™ Plant	Solution	Mx	227-201
GeneAll® GenEx™ Plant	Solution	Lx	227-301
GeneAll® GenEx™ Plant plus!	Solution	Sx	228-101
GeneAll® GenEx™ Plant plus!	Solution	Mx	228-250
GeneAll® GenEx™ Plant plus!	Solution	Lx	228-320


Sx On the basis of DNA purification from 100 mg plant tissue

Mx On the basis of DNA purification from 500 mg plant tissue

Lx On the basis of DNA purification from 2 g plant tissue

GeneAll® DirEx™

Procedures



Description

GeneAll® DirEx™ is for the preparation of template DNA for PCR in a single tube. It can be conveniently used for preparing the total DNA in just 8 minutes without the use of toxic chemicals such as phenol or chloroform. Template DNA can be prepared from various biological or forensic samples including whole blood, body fluids, tissues, cultured cells, tail snips, hair follicles, buccal swabs, cigarette butts and etc.

The procedure of DirEx™ consists of just two-steps, the lysis step and the inactivation step. The first is the lysis of samples with DirEx™ DNA preparation solution and proteinase K. And the following inactivation step inactivates the activity of proteinase K by heating. The prepared DNA solution can be used for the template of PCR or stored at -70 °C for long-term storage.


Features and Benefits

- Two steps, simple to use and user friendly kit for gDNA purification
- It takes 8 minutes to complete the whole procedure
- Optimized protocols for various samples such as cultured cell, animal tissue, hair, buccal swab, whole blood and dried blood spot
- Solution type
- No organic reagents, eco-friendly method
- Stable and consistent results
- Applicable to downstream such as PCR, Sequencing and RFLP etc.
- Proteinase K included

Component list

DirEx™ Solution
Proteinase K
PK Storage buffer
Manual

Stable and Consistent Result



PCR analysis of prepared DNA using **GeneAll® DirEx™** represents consistent result. Template DNA was isolated from blood (Lane 1~10).

Lane M : 100 bp ladder marker

Lane NC : Negative control

Single tube DNA preparation solution for PCR


DirEx™

Format : Solution

Sample size : 20 μ l whole blood
 2 x 10⁶ cultured cells
 10 mg animal tissue

Preparation time : 8 min

Stable and Consistent Result



PCR analysis was performed with prepared DNA using **GeneAll® DirEx™**.

Template DNA was isolated from CHO cells (Lane 1,2), RAW264.7 cells (Lane 3,4), Liver (Lane 5,6), Brain (Lane 7,8), Tail snips (Lane 9,10), Whole blood (Lane 11,12), Dried blood spot (Lane 13,14), Hair follicle (Lane 15,16), Hair shaft (Lane 17,18), Buccal swab (Lane 19,20), Cigarette butts (Lane 21~24).

Lane M : 100 bp ladder marker

Products	Type	Size	Cat. No.
GeneAll® DirEx™	Solution	50	250-050

For Genomic DNA Purification

Exgene™ / GenEx™ / DirEx™ Series

Exgene™ and GenEx™ series are designed for the purification of total DNA from a variety of sample sources. Exgene™ series provide fast and easy methods in convenient spin or vacuum column format and there are no need phenol extraction or alcohol precipitation. GenEx™ series provide convenient, scalable purification methods in the specially formulated buffer system. DirEx™ can be conveniently used for preparation of total DNA from various biological samples without the use of toxic chemical such as phenol or chloroform. Purified total DNA can be directly applicable in conventional PCR, real-time PCR, southern blotting, genotyping, RFLP and other downstream applications.

	Exgene™ Tissue SV (plus)*	Exgene™ Blood SV	Exgene™ Cell SV	Exgene™ Clinic SV	Exgene™ Genomic DNA micro	Exgene™ Plant SV	Exgene™ Soil DNA mini	GenEx™ Blood / Cell / Tissue**	GenEx™ Plant (plus)***	DirEx™
Sample Type	B C T									
Animal tissue	○		○	○	○			○		○
Body fluid		○	○	○	○			○	○	△
Bone					○					
Buccal swab	△	○	○	○	○			○	○	○
Buffy coat		○	○	○	△		△			
Callus						○			○	
Cultured cells	○	○	○	○	△			○	○	○
DNA virus		○	○	○	△			△		
Dried blood spot	△		○	○	○			△		
Fixed tissue	△		○	○	△			△		
Forensic sample					○					
Fungi						○	△		○	
Gram(-) bacteria	○		○	○	△			○	○	○
Gram(+) bacteria			○					△	△	△
Hair	△	○	○	○	○			△		○
Lichens							○			
Insect / worm	○		△	△	△			○		△
Mammalian whole blood	○*	○	○	○	○		○			○
Nail					○					○
Nucleated blood	△	○	○	○	△			△	△	△
Paraffin block	○		○	○	△			○		
Plant cells						○			○	
Plant tissue						○			○	
Rodent tails	○		○	○	△			○		○
Saliva		○	○	○	○			△		
Soil							○			
Sperm		○	○	○	○			△		
Urine	△		△	△	○			△		
Yeast			○					△	△	

○ Recommended / △ Suitable but not optimized and required additional protocol

* Exgene™ Tissue plus! provides the additional methods for the purification of total DNA from mammalian whole blood.

** GenEx™ series provide convenient, scalable purification methods in the specially formulated buffer systems.




*** GenEx™ Plant plus! kit has an additional feature, EzSep™ filter column for cleared supernatant

For Genomic DNA Purification

Exgene™ / GenEx™ / DirEx™ Series

	Exgene™ Tissue SV (plus!)	Exgene™ Blood SV	Exgene™ Cell SV	Exgene™ Clinic SV	Exgene™ Genomic DNA micro	Exgene™ Plant SV	Exgene™ Soil DNA mini	GenEx™ Blood / Cell / Tissue	GenEx™ Plant (plus!)	DirEx™
Specifications										
Format	Spin / Vacuum	Spin / Vacuum	Spin / Vacuum	Spin / Vacuum	Spin	Spin / Vacuum	Spin / Vacuum	Solution	Solution	Solution
Scale	mini Midi MAXI	mini Midi MAXI	mini MAXI	mini Midi MAXI	mini	mini Midi MAXI	mini	Sx Lx	Sx Mx Lx	50 ml
Starting sample	25 mg 100 mg 250 mg	300 µl 2 ml 10 ml	5 x 10 ⁶ cells 1 x 10 ⁸ cells	25 mg 100 mg 250 mg	#	100 mg 400 mg 1000 mg	500 mg	#	#	#
Typical yield	5 ~ 50 µg 20 ~ 150 µg 80 ~ 400 µg	4 ~ 20 µg 20 ~ 80 µg 80 ~ 400 µg	10 ~ 50 µg 80 ~ 500 µg	5 ~ 50 µg 20 ~ 150 µg 80 ~ 400 µg	#	4 ~ 40 µg 10 ~ 150 µg 40 ~ 300 µg	#	#	#	#
Preparation time	25 ~ 220 min 40 ~ 250 min 40 ~ 250 min	20 ~ 30 min 40 ~ 55 min 40 ~ 55 min	30 ~ 120 min 60 ~ 240 min	25 ~ 220 min 40 ~ 250 min 40 ~ 250 min	25 ~ 220 min	< 40 min < 1 hour < 1 hour	25 min	25 ~ 90 min 25 ~ 90 min	< 30 min < 60 min < 70 min	< 8 min

Typical yield depends on the type and size of sample.



Elution Tip!

The yield and concentration of purified DNA depending on the elution method. DNA was prepared from 200 µl of bovine whole blood. Each preparation was exactly identical except the elution method; Elution was performed 3 times per column with 100 µl (Lane 1 ~ 3) and 200 µl (Lane 5 ~ 7) and 2 times per column with 300 µl (Lane 9 ~ 10) of fresh buffer AE. At the same time, another elution was carried out 3 times (100 µl and 200 µl) and 2 times (300 µl) by recursive use of the eluate instead of fresh buffer AE (Lane 4, 8, 11). Total 11 eluates purified from 6 samples were resolved on 0.8 % agarose gel to visualize (upper left) and its concentration (lower left) and total yield (lower right) was measured by spectrophotometric analysis.

Coming Soon.....

Exgene™ Stool DNA mini

For more information about products, visit www.geneall.com


4. RNA Purification System



Hybrid-R™	50
Hybrid-R™ Blood RNA	52
Hybrid-R™ miRNA	54
RiboEx™	56
RiboEx™ LS	58
Ribospin™	60
Ribospin™ vRD (plus!)	62
Ribospin™ Plant	64
Riboclear™ (plus!)	66
Allspin™	68
Selection Guide for RNA Purification	70

GeneAll[®] Hybrid-R[™]

Procedures



Component list

RiboEx[™]

GeneAll[®] column type B

2 ml collection tube

1.5 ml collection tube

Buffer RB1

Buffer SW1

Buffer RNW

RNase-free water

Protocol Handbook

Description


GeneAll[®] Hybrid-R[™] provides an easy and rapid method for the isolation of highly purified total RNA from samples of human, animal, plant, yeast, bacterial and viral origin. Hybrid-R[™] eliminates alcohol precipitation by binding of RNA with column, allowing rapid and convenient preparation from a large number of samples simultaneously.

Hybrid-R[™] can yield up to 500 µg depending on the type of tissue sample used and complete all process to prepare total RNA in just 30 minutes. The purified total RNA is suitable for the isolation of mRNA, northern blotting, dot blotting, in vitro translation, cloning, RT-PCR, RNase protection assays and other analytical procedures.

Features and Benefits

- Preparation time : ~ 30 minutes
- Accurate and consistent yield from animal tissue, cultured cell line, plant, E. coli and various biological samples
- High purity and yield
- No Genomic DNA contamination
- No ethanol precipitation
- Ready for use in RT-PCR, northern blotting, dot blotting, in vitro translation, molecular cloning, real-time PCR, RNase protection assays and other analytical procedures

RNA Purification Result



Total RNA was purified from E. coli ($OD_{600} \approx 1.8$) using several RNA extraction kits of different companies. E. coli cells were taken to the total RNA purification.


The purified total RNA was loaded on a 1 % formaldehyde gel.

M : 0.5 ~ 10 kb RNA ladder

Lane 1 : Total RNA from Hybrid-R[™]

Lane 2 : Total RNA from supplier A

Lane 3 : Total RNA from supplier B



Total RNA was purified from rat liver using several RNA extraction kits of different companies. 100 mg / 1.2 ml was taken to the total RNA purification.

The purified total RNA was loaded on a 1 % formaldehyde gel.

M : 0.5 ~ 10 kb RNA ladder

Lane 1 : Total RNA from Hybrid-R[™]

Lane 2 : Total RNA from supplier A

Lane 3 : Total RNA from supplier B

Lane 4 : Total RNA from supplier C

For the isolation of total RNA from tissues and cultured cells

Hybrid-R™



Format : Spin column type B
with 1.5 ml / 2 ml collection tubes


Sample size : ~ 100 mg tissue or 1×10^7 cells

Application volume : ~ 700 μl

Minimum elution volume : 30 μl

Binding capacity : ~ 500 μg

RT-PCR Result



Total RNA was purified from various samples using Hybrid-R™.
And then, the cDNA was synthesized by reverse transcriptase.
The cDNA was amplified by PCR and confirmed by electrophoresis.

M : Lambda-HindIII

Lane 1 : PCR of E. coli cDNA

Lane 2 : PCR of Rat kidney cDNA


Lane 3 : PCR of Rat liver cDNA

Lane 4 : PCR of Rat heart cDNA

Products	Type	Size	Cat. No.
GeneAll® Hybrid-R™	Spin	100	305-101

GeneAll[®] Hybrid-R[™] Blood RNA

Procedures



Component list

GeneAll[®] EzPure[™] filter
 GeneAll[®] Column type W
 2 ml collection tube
 1.5 ml collection tube
 RiboEx[™] LS
 Buffer RB1
 Buffer RBW
 Buffer RNW
 RNase-free water
 Protocol Handbook

Description

GeneAll[®] Hybrid-R[™] Blood RNA is a complete kit with ready-to-use reagent for the isolation of total RNA from up to 0.25 ml whole blood sample.

This kit utilizes the lysis method of RiboEx[™] LS which has a powerful ability of cell-lysis and the purification method based on glassfiber membrane technology. Fast and convenient procedure of Hybrid-R[™] Blood RNA takes only 30 minutes for complete preparation of pure RNA.

Whole blood sample is homogenized and lysed in RiboEx[™] LS, a monophasic solution containing phenol and guanidium salt, which rapidly lyse cells and inactivates nucleases. In conventional methods, the erythrocytes of mammalian blood which does not contain nuclei (and therefore, RNA either) should be removed by pre-treatment such as osmotic lysis for the separation of leukocytes from whole blood. This additional treatments increase the experiment time and the possibility of RNA-breakage, followed by decline of RNA-quality.

Hybrid-R[™] Blood RNA does not need the additional treatment of blood sample, and whole blood is lysed in RiboEx[™] LS in just one step.

Then addition of chloroform brings about a separation of the lysate into aqueous and organic phases. After phase-separating, DNA and protein remains in the interphase and the organic phase respectively but released RNA exists in the aqueous phase.

The aqueous phase is picked and applied to a EzPure[™] filter to eliminate small amount of contaminated DNA and other blood contaminants. The passed-through is mixed with buffer RB1, RNA binding buffer and then the mixture is applied to a mini spin column. After a series of washing with buffer RBW and RNW, pure RNA can be eluted by RNase-free water.

Hybrid-R[™] Blood RNA is suitable for RNA preparation from 0.1 ml to 0.25 ml mammalian whole blood. The typical yield is 3 µg per 0.25 ml whole blood. The purified RNA can be applicable for the isolation of Poly A⁺ RNA, northern blotting, dot blotting, in vitro translation, cloning, RT-PCR, RPA and other analytical procedures.

Features and Benefits

- Preparation time : ~ 30 minutes
- Accurate and consistent yield from whole blood
- High purity and yield
- Sample size : 100 ~ 250 µl / prep
- Instant use : No need of additional materials
- No ethanol precipitation
- No Genomic DNA contamination

For the isolation of total RNA from whole blood

Hybrid-R™ Blood RNA



Format : Spin column type W
with 2 ml collection tubes


Sample size : 100 ~ 250 μl whole blood

Application volume : ~ 700 μl

Minimum elution volume : 30 μl

Binding capacity : 100 μg

Comparison Data



Total RNA was extracted from whole blood using several RNA extraction kits of different companies. The extracted total RNA was loaded on a 1 % formaldehyde gel.


Lane 1 : Total RNA from Hybrid-R™ Blood RNA for 250 μl of whole blood

Lane 2 : Total RNA from supplier A for 500 μl of whole blood

Lane 3 : Total RNA from supplier B for 500 μl of whole blood

Lane 4 : Total RNA from supplier C for 250 μl of whole blood

Genomic DNA contamination and RT-PCR result



As analysis of genomic DNA contamination, PCR for human beta-actin was performed with eluates, purified from whole blood using several kits of other companies.


M : 1 Kb ladder

Lane 1, 2 : PCR of the eluate from Hybrid-R™ Blood RNA

Lane 3, 4 : PCR of the eluate from supplier A

Lane 5, 6 : PCR of the eluate from supplier B

Lane 7, 8 : PCR of the eluate from supplier C



Total RNA was extracted from whole blood using Hybrid-R™ Blood RNA and other supplier kits. And then the cDNA was synthesized by reverse transcriptase. The cDNA was amplified by human beta-actin primer and confirmed by electrophoresis.

M : 1 Kb ladder

Lane 1 : PCR of cDNA from Hybrid-R™ Blood RNA


Lane 2 : PCR of cDNA from supplier A

Lane 3 : PCR of cDNA from supplier B

Products	Type	Size	Cat. No.
GeneAll® Hybrid-R™ Blood RNA	spin	50	315-150

GeneAll[®] Hybrid-R[™] miRNA

Procedures



Component list

GeneAll[®] Column type B (red ring)
(with collection tube)
GeneAll[®] Column type W (blue ring)
(with collection tube)
2 ml collection tube
1.5 ml collection tube
RiboEx[™]
Buffer SW1
Buffer RBW
Buffer RNW
RNase-free water
Protocol Handbook

Description


In recent years, interest in small RNA, such as siRNA and miRNA which are related to research of gene regulation, has expanded. There are many commercial kits for total RNA preparation, but most of these are focused on preparation of large RNA longer than 200 nucleotides. Because both siRNA and miRNA are between 15 ~ 30 nucleotides in length, the need of specially optimized kit for small RNA (< 200 nucleotides) is growing rapidly. Hybrid-R[™] miRNA is designed for purification of large and small RNA separately from culture cells or animal tissues and co-purification in a single tube is also available by modified protocol. This kit utilizes the lysis method of RiboEx[™] which has a powerful ability of lysis and the purification method based on glassfiber membrane technology.

Samples are homogenized in RiboEx[™], a monophasic solution containing phenol and guanidium salt, which rapidly lyse cells and inactivates nucleases. Addition of chloroform brings about a separation of the lysate into aqueous and organic phases. Total RNA locates in the aqueous phase while DNA and protein remain in the interphase and organic phase. Large and small RNA in the aqueous phase is selectively bound to column type B and type W respectively. The column type B selectively adsorbs the RNA larger than 200 nucleotides in length, while the column type W specifically holds the RNA smaller than 200 nucleotides in length.

To purify large RNA, the aqueous phase is mixed with ethanol and the mixture is applied to a column type B. After centrifugation, large RNA is bound to membrane and the mixture containing small RNA goes into collection tube through the membrane. The membrane is washed away by two wash buffer (SW1 and RNW) and purified large RNA is eluted from the membrane by RNase-free water.

To purify small RNA, the pass-through come from the binding of large RNA is mixed with ethanol and then applied to a column type W. After washing with buffer RBW and RNW, small RNA is eluted by RNase-free water. The procedure of Hybrid-R[™] miRNA takes only 30 minutes for complete preparations of pure RNA. The purified RNA is suitable for the isolation of Poly A⁺ RNA, northern blotting, dot blotting, in vitro translation, cloning, RT-PCR, RPA and other analytical procedures.

Real-time PCR Result of miR-24 from small RNA




Real-time PCR was performed with purified small (micro) RNA using Hybrid-R[™] miRNA kit. Small RNA was extracted from CHO cell, RAW264.7 cell and rat heart and liver. And then RT-PCR of miR-24 was performed using miScript PCR system (Qiagen). Amplified miR-24 was detected by 7500 Real-Time PCR system (Applied Biosystems).

For purification of large and small RNA separately from cultured cells or animal tissues

Features and Benefits

- Preparation time : ~ 30 minutes
- Stable and consistent yield
- High purity and yield
- Perfect separation of small RNA fragment
- Sample size : Up to 50 mg tissue or up to 1×10^7 cultured cells
- Recovery range : Large RNA : > 200 nucleotides
Small RNA : < 200 nucleotides
- Instant use : No need of additional materials
- No ethanol precipitation
- No Genomic DNA contamination
- Ready for use in northern blotting, dot blotting, in vitro translation, cloning, RT-PCR, RPA and other analytical procedures

Experimental Result



Large and small RNA was extracted from CHO (chinese hamster ovary) cell, RAW264.7 cell and rat lung tissue using Hybrid-R™ miRNA. The purified large RNA was loaded on a 1% formaldehyde gel and small RNA was loaded on a 1% agarose gel.

M1 : 0.5 ~ 10 kb RNA ladder

M2 : Lambda-HindIII

Lane 1 : Large RNA from CHO cell


Lane 4 : Small RNA from CHO cell

Lane 2 : Large RNA from RAW264.7 cell

Lane 5 : Small RNA from RAW264.7 cell

Lane 3 : Large RNA from rat lung

Lane 6 : Small RNA from rat lung



miRNA was extracted using several miRNA extraction kits of different companies. The extracted miRNA was loaded on a 15 % urea-acrylamide gel.

Lane 1 : miRNA from Hybrid-R™ miRNA for CHO cell

Lane 2 : miRNA from Hybrid-R™ miRNA for RAW264.7 cell


Lane 3 : miRNA from Hybrid-R™ miRNA for rat lung

Lane 4 : miRNA from supplier A for CHO cell

Lane 5 : miRNA from supplier A for RAW264.7

Lane 6 : miRNA from supplier A for rat lung

RT-PCR result of large RNA



Large RNA was purified from CHO cell, RAW264.7 cell and rat lung tissue using Hybrid-R™ miRNA.

And then cDNA was synthesized by reverse transcriptase. The cDNA was amplified by beta-actin primer and confirmed by electrophoresis.

M : 1kb ladder

Lane 1 : PCR of cDNA from CHO cell


Lane 2 : PCR of cDNA from RAW264.7

Lane 3 : PCR of cDNA from rat lung

Products	Type	Size	Cat. No.
GeneAll® Hybrid-R™ miRNA	spin	50	325-150

GeneAll[®] RiboEx[™]

Procedures



Component list

GeneAll[®] RiboEx[™] Total RNA
Isolation Solution
Manual


Description

GeneAll[®] RiboEx[™] is a complete kit with ready-to-use reagents for the isolation of total RNA from samples of human, animal, plant, yeast, bacterial and viral origin. RiboEx[™] is based on the disruption of cells in guanidine salt / detergent solution, followed by organic extraction and alcohol precipitation of the RNA, and it allows simultaneous processing of a large number of samples. RiboEx[™] can yield up to 10 μg / 1 mg tissue or up to 30 μg / 1×10^6 cultured cells of highly purified total RNA. The resulting total RNA is suitable for the isolation of poly A⁺ RNA, northern blotting, dot blotting, in vitro translation, cloning, RT-PCR, RNase protection assays and other analytical procedures.

Features and Benefits


- Format : Monophase solution type
- Sample size : Up to 100 mg tissue / 1 ml RiboEx[™]
Up to 1×10^7 cells / 1 ml RiboEx[™]
- Preparation time : 50 ~ 65 minutes
- Typical yield : Up to 10 μg / 1 mg tissue
Up to 30 μg / 1×10^6 cultured cells
- High purity : $\text{OD}_{260/230} > 2.0$, $\text{OD}_{260/280} > 1.8$
- Accurate and consistent yield from animal tissue, cultured cell line, plant, E. coli and various biological samples
- Accurate and easy phase separation
- Ready for use in RT-PCR, northern blotting, dot blotting, in vitro translation, molecular cloning, real-time PCR, RNase protection assays and other analytical procedures

RNA Purification Result



Total RNA was purified from rat brain using several RNA extraction kits of different companies. 100 mg / 1.2 ml was taken to the total RNA purification. The purified total RNA was loaded on a 1 % formaldehyde gel.

Lane 1 : Total RNA from RiboEx[™]
Lane 2 : Total RNA from Supplier A
Lane 3 : Total RNA from Supplier B
Lane 4 : Total RNA from Supplier C



Total RNA was purified from E. coli DH5α using several RNA extraction kits of different companies. E. coli cells were taken to the total RNA purification. The purified total RNA was loaded on a 1 % formaldehyde gel.


Lane 1 : Total RNA from RiboEx[™]
Lane 2 : Total RNA from Supplier A
Lane 3 : Total RNA from Supplier B
Lane 4 : Total RNA from Supplier C

For total RNA isolation from various samples

Total RNA yield from various starting materials using RiboEx™

Material	Sample type	Amount	Yield of RNA
Cell Lines	RAW 264.7	1 x 10 ⁶ cells	~ 28 µg
Animal Tissue	Liver	1 mg	~ 10 µg
	Spleen	1 mg	~ 10 µg
	Kidney	1 mg	~ 4 µg
	Brain	1 mg	~ 1.5 µg
Gram(-) Bacteria	E. coil	O.D ₆₀₀ = 1.8 (1.5 ml pellet)	~ 60 µg

RT-PCR Result



Total RNA was purified from Mouse ES cell using RiboEx™ and supplier A kits. And then, the cDNA was synthesized by reverse transcriptase. The cDNA was amplified by PCR and confirmed by electrophoresis.

M : 1 Kb ladder

Lane 1, 3, 5 : PCR of cDNA from RiboEx™

Lane 2, 4, 6 : PCR of cDNA from supplier A

Lane 1, 2 : amplified by β-actin primer

Lane 3, 4, 5, 6 : amplified by Oct 4 primer

Fig I, II Total RNA was purified from 293 cell using RiboEx™ and supplier A kits. And then the cDNA was synthesized by reverse transcriptase. The cDNA was amplified by PCR and confirmed by electrophoresis.


Lane A : supplier A kit

Lane B : RiboEx™

Fig I : amplified by GAPDH primer

Fig II : amplified by Hif-1 primer

Real-Time PCR amplification




Real-time PCR was performed with purified total RNA using RiboEx™ total RNA isolation kit. Total RNA was extracted from rat atrium. And then the cDNA was synthesized through reverse transcriptase. Reference was confirmed by GAPDH primer and target gene was confirmed by ET-1 primer in the experimental and control group.

Products	Type	Size	Cat. No.
GeneAll® RiboEx™	Solution	100	301-001
GeneAll® RiboEx™	Solution	200	301-002

GeneAll[®] RiboEx[™] LS

Procedures



Component list

GeneAll[®] RiboEx[™] LS Total RNA Isolation Solution
Manual

Description

GeneAll[®] RiboEx[™] LS is a complete kit with ready-to-use reagents for the isolation of total RNA from various liquid samples. RiboEx[™] LS is a concentrated form of RiboEx[™] and this allows that liquid samples can be processed more easily with it, while RiboEx[™] is more suitable for solid samples and pelleted cells. RiboEx[™] LS is a mono-phasic solution containing phenol and guanidine salt, which rapidly lyse cells and inactivates nucleases. Addition of chloroform brings about a separation of the homogenate in aqueous and organic phases. RNA locates in the aqueous phase while DNA and protein remain in the interphase and organic phases. The aqueous phase including RNA is mixed with isopropanol and the RNA which is precipitated by centrifuging. The purified total RNA is suitable for RT-PCR, northern blotting, dot blotting, in vitro translation, molecular cloning, real-time PCR, RNase protection assays and other analytical procedures.

Features and Benefits

- Format : Monophase solution type
- Sample size : Up to 0.25 ml liquid sample / 0.75 ml RiboEx[™] LS
Up to 100 mg tissue / 0.75 ml RiboEx[™] LS
- Preparation time : 50 ~ 65 minutes
- Typical yield : Up to 30 μ g / 1×10^6 cultured cells
Up to 10 μ g / 1 mg tissue
- High purity : $A_{260}/A_{230} > 2.0$, $A_{260}/A_{280} > 1.8$
- Accurate and consistence yield
- Ready for use in RT-PCR, northern blotting, dot blotting, in vitro translation, molecular cloning, real-time PCR, RNase protection assays and other analytical procedures

Genomic DNA Contamination Test and RT-PCR Result




Fig. I Genomic DNA contamination was tested by PCR. Eluate, including total RNA of RAW264.7 cell, from several RNA extraction kits of different companies was the template of PCR and amplified by beta-actin primer.

M : Lambda-HindIII

Lane 1 : PCR of the eluate from supplier A for liquid sample

Lane 2 : PCR of the eluate from RiboEx[™] LS

Lane 3 : PCR of the eluate from supplier B for liquid sample

Fig. II Total RNA was extracted from RAW264.7 cell using RiboEx[™] LS and supplier kits. And then, the cDNA was synthesized by reverse transcriptase. The cDNA was amplified by PCR and confirmed by electrophoresis.

Lane 4 : PCR of cDNA from supplier A for liquid sample

Lane 5 : PCR of cDNA from RiboEx[™] LS


Lane 6 : PCR of cDNA from supplier B for liquid sample

For total RNA isolation from various liquid samples

Total RNA yield from various starting materials using RiboEx™ LS

Material	Sample type	Amount	Yield of RNA
Cell Lines	RAW 264.7	1 x 10 ⁶ cells	~ 28 µg
Animal Tissue	Liver	1 mg	~ 10 µg
	Spleen	1 mg	~ 10 µg
	Kidney	1 mg	~ 4 µg
	Brain	1 mg	~ 1.5 µg
Blood	Whole human or animal blood	1 mg	~ 1.5 µg
Gram(-) Bacteria	E. coil	O.D ₆₀₀ = 1.8 (1.5 ml pellet)	~ 60 µg

Comparison Data




Total RNA was extracted from E. coli DH5α using several RNA extraction kits of different companies. The extracted total RNA was loaded on a 1 % formaldehyde gel.

Lane 1 : Total RNA from Supplier A for liquid sample

Lane 2 : Total RNA from RiboEx™ LS

Lane 3 : Total RNA from Supplier B for liquid sample




Total RNA was extracted from CHO (chinese hamster ovary) cell using several RNA extraction kits of different companies. The extracted total RNA was loaded on a 1 % formaldehyde gel.

Lane 1 : Total RNA from Supplier A for liquid sample

Lane 2 : Total RNA from RiboEx™ LS

Lane 3 : Total RNA from Supplier B for liquid sample




Total RNA was extracted from heart tissue of rat using several RNA extraction kits of different companies. The extracted total RNA was loaded on a 1 % formaldehyde gel.

Lane 1 : Total RNA from Supplier A for liquid sample

Lane 2 : Total RNA from RiboEx™ LS

Lane 3 : Total RNA from Supplier B for liquid sample



Total RNA was extracted from whole blood of rat using several RNA extraction kits of different companies. The extracted total RNA was loaded on a 1 % formaldehyde gel.

Lane 1 : Total RNA from Supplier A for liquid sample


Lane 2 : Total RNA from RiboEx™ LS

Lane 3 : Total RNA from Supplier B for liquid sample

Products	Type	Size	Cat. No.
GeneAll® RiboEx™ LS	Solution	100	302-001
GeneAll® RiboEx™ LS	Solution	200	302-002

GeneAll[®] Ribospin[™]

Procedures



Component list

- GeneAll[®] column type F
- 2 ml collection tube
- 1.5 ml collection tube
- Buffer LYS
- Buffer GW1
- Buffer RNW
- RNase-free water
- Manual

Description

Ribospin[™] provides a convenient method for isolation of total RNA from cell and tissue samples. Ribospin[™] procedures employed the glassfiber membrane technology for the fastest and the most convenient of high purity RNA isolation, instead of conventional alcohol precipitation or phenol / chloroform extraction.

Whole procedure takes only 15 minutes and the eluates are suitable for RT-PCR or any downstream application without further manipulation.

Features and Benefits

- Glassfiber membrane technology
- Sample size : Up to 25 mg tissue or up to 5 x 10⁶ cultured cells
- Typical yield : Up to 20 µg / 1 x 10⁶ cultured cells
Up to 60 µg / 10 mg liver tissue
- High purity : A₂₆₀/A₂₃₀ > 2.0, A₂₆₀/A₂₈₀ > 1.8
- Preparation time : ~ 15 minutes
- Stable and consistent yield
- No phenol / chloroform extraction
- No ethanol precipitation
- Ready for use in RT-PCR, northern blotting, dot blotting, in vitro translation, molecular cloning, real-time PCR, RNase protection assays and other analytical procedures

The yield of total RNA may vary depending on the tissue or cells from which it is obtained.

Material	Sample type	Amount of starting material	Yield of Total RNA
Cultured cell	CHO	1 x 10 ⁶ cells	~ 15 µg
	RAW 264.7	1 x 10 ⁶ cells	~ 20 µg
Tissue	Liver	10 mg	~ 60 µg
	Kidney	10 mg	~ 30 µg
	Spleen	10 mg	~ 35 µg
E. coli	DH5α	O.D 1.5 (2 ml pellet)	~ 10 µg

For total RNA isolation from animal tissues and cultured cells

Ribospin™



Format : Spin column type F
with 1.5 ml / 2 ml collection tubes


Sample size : ~ 25 mg or 5×10^6 cells

Application volume : ~ 750 μl

Minimum elution volume : ~ 40 μl

Binding capacity : ~ 100 μg

Downstream Application Test




Total RNA was extracted from RAW264.7 cell using several RNA extraction kits of different companies. The extracted RNA was loaded on a 1% formaldehyde gel.

M : 0.5 ~ 10 kb RNA ladder

Lane 1 : Total RNA from supplier A

Lane 2 : Total RNA from supplier B

Lane 3 : Total RNA from Ribospin™




Total RNA was extracted from CHO (chinese hamster ovary) cell using several RNA extraction kits of different companies. The extracted RNA was loaded on a 1% formaldehyde gel.

M : 0.5 ~ 10 kb RNA ladder

Lane 1 : Total RNA from supplier A

Lane 2 : Total RNA from supplier B

Lane 3 : Total RNA from Ribospin™




Total RNA was extracted from RAW264.7 cell using Ribospin™ and supplier kits. And then the cDNA was synthesized by reverse transcriptase. The cDNA was amplified by PCR and confirmed by electrophoresis.

M : Lambda-HindIII

Lane 1 : PCR of cDNA from supplier A

Lane 2 : PCR of cDNA from supplier B

Lane 3 : PCR of cDNA from Ribospin™



Total RNA was extracted from liver tissue of rat using Ribospin™ and supplier kits. And then the cDNA was synthesized by reverse transcriptase. The cDNA was amplified by PCR and confirmed by electrophoresis.

M : Lambda-HindIII

Lane 1 : PCR of cDNA from supplier A


Lane 2 : PCR of cDNA from supplier B

Lane 3 : PCR of cDNA from Ribospin™

Products	Type	Size	Cat. No.
GeneAll® Ribospin™	Spin	50	304-150

GeneAll[®] Ribospin[™] vRD (plus!)

Procedures



Description

Ribospin[™] vRD provides a convenient method for isolation of RNA and DNA from cell-free fluid, cell-culture supernatant, plasma, serum, swab, urine, and virus-infected samples. Ribospin[™] vRD procedures employed the glassfiber membrane technology for the fastest and the most convenient of high purity RNA and DNA isolation, instead of conventional alcohol precipitation or phenol/chloroform extraction. Ribospin[™] vRD buffer system provides the effective binding condition of RNA and DNA to glassfiber membrane through mix with lysis and binding buffers. And then the impurities on the membrane are washed away by two different wash buffers. At last, pure RNA and DNA are eluted by nuclease-free water. Whole procedure may take only 20 minutes and the eluate is suitable for PCR, RT-PCR or any downstream application without further manipulation.

Ribospin[™] vRD plus! kit offers carrier RNA for purification of nucleic acid from very small amounts of sample.

Features and Benefits

- Glassfiber membrane technology
- Sample size : up to 300 μl
- Preparation time : ~ 20 minutes
- Stable and consistent yield
- No phenol / chloroform extraction
- No ethanol precipitation
- Ready for use in PCR, RT-PCR, real-time PCR and other analytical procedures

Component list

GeneAll[®] column type V
 2 ml collection tube
 1.5 ml collection tube
 Buffer VL
 Buffer RB1
 Buffer RBW
 Buffer RNW
 Nuclease-free water
 Carrier RNA (plus! only)
 Manual

Ribospin[™] vRD (plus!)



Format : Spin column type V
 with 1.5 ml / 2 ml collection tubes

Sample size : ~ 300 μl

Application volume : ~ 800 μl

Minimum Elution volume : ~ 30 μl

Binding capacity : ~ 100 μg

For viral RNA / DNA isolation from various samples

Experimental Result

* Amplification test of HPIV (human parainfluenza virus) RNA



Fig I : PCR result using HPIV 1 specific primer



Fig II : PCR result using HPIV 2 specific primer




Fig III : PCR result using HPIV 3 specific primer

Viral RNA was purified from HPIV (human parainfluenza virus) 1, 2, 3 using Ribospin™ vRD.
And then the cDNA was synthesized by reverse transcriptase. The cDNA was amplified by PCR and confirmed by electrophoresis.

M : 100 bp ladder

Lane 1 ~ 3 : First PCR result

Lane 4 ~ 6 : Nest PCR result

Lane 1, 4 : HPIV 1

Lane 2, 5 : HPIV 2

Lane 3, 6 : HPIV 3

* Amplification test of HSV-1 (Herpes simplex virus) DNA




Total nucleic acid was extracted from cells infected by HSV-1 (DNA virus) and HSV-1 samples using Ribospin™ vRD. The DNA of HSV-1 was amplified by PCR and confirmed by electrophoresis.

M : Lambda-HindIII

Lane 1 : PCR of DNA from infected cell

Lane 2, 3 : PCR of DNA from HSV-1 sample

Lane 4 : Negative control



Total DNA was extracted from gradually diluted HSV-1 sample using Ribospin™ vRD. And then the DNA of HSV-1 was amplified by PCR and confirmed by electrophoresis.

M : Lambda-HindIII

Lane 1 : PCR of DNA extracted from 6×10^4 pfu HSV-1

Lane 2 : PCR of DNA extracted from 6×10^3 pfu HSV-1

Lane 3 : PCR of DNA extracted from 6×10^2 pfu HSV-1

Lane 4 : PCR of DNA extracted from 6×10 pfu HSV-1

Lane 5 : PCR of DNA extracted from 6 pfu HSV-1


Lane 6 : Negative control of a purification procedure

Lane 7 : Negative control

Products	Type	Size	Cat. No.
GeneAll® Ribospin™ vRD	Spin	50	302-150
GeneAll® Ribospin™ vRD plus!	Spin	50	312-150

GeneAll[®] Ribospin[™] Plant

Procedures



Component list

GeneAll[®] EzPure[™] filter (yellow)
(with collection tube)
GeneAll[®] column type W (blue ring)
(with collection tube)
1.5 ml collection tube
DNase I
Buffer RPL
Buffer REL
Buffer RBW
Buffer RNW
RNase-free water
Buffer DRB
Protocol Handbook

Description

Ribospin[™] Plant is specially designed for purification of total RNA from various plant tissues such as leaves, stems, roots and picky plant samples. This kit provides the optimized buffer and spin column, which is effective in removing polysaccharides and polyphenolic compounds and isolating intact plant RNA. All components of Ribospin[™] Plant are ready to use, so any further preparation for experiment is not required.

The procedure of Ribospin[™] Plant begins with the disruption of sample in liquid nitrogen using mortar and pestle. The disrupted sample can be lysed in buffer RPL or REL. In most case, buffer RPL is the best buffer for lysis. However in some plant samples, solidification of lysate can be occurred with buffer RPL due to endosperm of seed or peculiar metabolites and this can be avoided by using buffer REL as alternative for buffer RPL.

Most impurities except RNA in the lysate are eliminated by filtration through EzPure[™] filter and then the passed-through lysate is mixed with ethanol to adjust binding condition. Total RNA including a little impurity is bound to the membrane of spin column type W while the mixture is passing through. Survived genomic DNA can be exterminated by on-column DNase I treatment at this step. After a series of washing step using buffer RBW and RNW, plant total RNA is eluted by RNase-free water. Whole procedure of Ribospin[™] Plant takes only 25 minutes. The purified RNA is suitable for cDNA synthesis, RT-PCR, northern blotting and other analytical procedure.

Features and Benefits

- Glassfiber membrane technology
- Sample size : Up to 100 mg plant tissue
- Including DNase I and treatment step
- Typical yield of RNA : ~ 50 µg / 100 mg tissue
- High purity : A_{260}/A_{280} 1.8~2.2, A_{260}/A_{230} >2.0
- Preparation time : ~ 25 minutes
- No phenol / chloroform extraction
- No ethanol precipitation
- Ready for use in RT-PCR, northern blotting, dot blotting, in vitro translation, molecular cloning, real-time PCR, RNase protection assay and other analytical procedures

For total RNA isolation from various plant samples

Ribospin™ Plant



Format : Spin column type W with 1.5 ml / 2 ml collection tubes,
EzPure™ filter with 2 ml collection tubes

Sample size : ~100 mg plant tissue

Maximum loading volume of


EzPure™ filter : ~ 600 µℓ

Application volume of spin column : ~ 700 µℓ

Minimum elution volume : 30 µℓ


Binding capacity : ~100 µg

RNA Purification Result



Total RNA was extracted from a wide variety of plant species using Ribospin™ Plant. The extracted RNA was loaded on a 1% formaldehyde gel.

Lane 1 : Leaf RNA from *Pinus densiflora*
Lane 2 : Leaf RNA from *Crassula ovata*
Lane 3 : Leaf RNA from *Citrus grandis* Osbek
Lane 4 : Leaf RNA from *Diospyros kaki*
Lane 5 : Leaf RNA from *Zea mays*
Lane 6 : Leaf RNA from *Lycopersicum esculentum*
Lane 7 : Leaf RNA from *Nicotiana tabacum*
Lane 8 : Leaf RNA from *Lactuca sativa*
Lane 9 : Leaf RNA from *Cucumis sativus* L
Lane 10 : Root RNA from *Plantago asiatica*
Lane 11 : Root RNA from *Nicotiana tabacum*
Lane 12 : Fruit RNA from *Citrus grandis* Osbek



Total RNA was purified from *Pinus densiflora* by Ribospin™ Plant. And the cDNA was synthesized by reverse transcriptase.


The cDNA was amplified by PCR and confirmed on a 1% agarose gel containing ethidium bromide.

Lane 1 : cDNA 1 µℓ
Lane 2 : cDNA 2 µℓ
Lane 3 : cDNA 3 µℓ
Lane 4 : cDNA 4 µℓ

Products	Type	Size	Cat. No.
GeneAll® Ribospin™ Plant	spin	50	307-150

GeneAll[®] Riboclear[™] (plus!)

Procedures



Description

Riboclear[™] provides an easy and rapid method for RNA cleanup or concentration from various RNA samples in just 6 minutes.

Riboclear[™] eliminates alcohol precipitation by binding of RNA with column, allowing rapid and convenient preparation from various samples simultaneously.

Purified RNA with Riboclear[™] series are free of salts and enzymes in yields reaching 95 % and are suitable for dot blotting, in vitro translation, cloning, RT-PCR, RNase protection assays and other analytical procedures. Riboclear[™] plus! kit provides DNase for removal of DNA and micro column for concentration of total RNA


Features and Benefits

- Preparation time : ~ 6 minutes
- Stable and consistent yield
- High recovery rate : ~ 95 %
- Instant use : No need of additional materials
- No use of organic solvents
- No ethanol precipitation
- Complete removal of salts and enzymes
- Ready for use in RT-PCR, northern blotting, dot blotting, in vitro translation, molecular cloning, real-time PCR, RNase protection assays and other analytical procedures

Component list

GeneAll[®] column type W
 2 ml collection tube
 1.5 ml collection tube
 Buffer MS
 Buffer RNW
 RNase-free water
GeneAll[®]
 Micro column type S (plus! only)
 DNase I (plus! only)
 Manual

Consistency Test



The consistency of the purified RNA using Riboclear[™] was confirmed by electrophoresis.

M : Lambda-HindIII

Lane 1 : Total RNA from Hybrid-R[™]

Lane 2 ~ 11 : The purified RNA from Riboclear[™]

For RNA cleanup from various RNA samples

Riboclear™ (plus!)



Format : Spin column type W (Micro column type S)
with 1.5 ml / 2 ml collection tubes

Sample size : ~ 100 $\mu\ell$

Recovery Rate : ~ 95%


Preparation time : ~ 6 minutes

Application volume : ~ 800 $\mu\ell$

Minimum Elution volume : 30 $\mu\ell$

Binding capacity : ~ 100 μg

Downstream Application Test



The purified RNA using Riboclear™ plus!
And then the cDNA was synthesized by reverse transcriptase.
The cDNA was amplified by PCR and confirmed by electrophoresis.
M : Lambda-HindIII

Lane 1 : PCR of undigested RNA eluate

Lane 2 ~ 4 : PCR of RNA eluate digested by DNase I


Lane 5 : RT-PCR of the product of Lane 1

Lane 6 ~ 8 : RT-PCR of the product of Lane 2, 3 and 4

Products	Type	Size	Cat. No.
GeneAll® Riboclear™	Spin	50	303-150
GeneAll® Riboclear™ plus!	Spin	50	313-150

GeneAll® Allspin™

Procedures



Component list

- GeneAll® column type B (red ring) (with collection tube)
- GeneAll® column type W (blue ring) (with collection tube)
- 2 ml collection tube without column
- 1.5 ml microcentrifuge tube
- Buffer CTL
- Buffer GW1
- Buffer BW
- Buffer RNW
- Buffer AE
- RNase-free water
- Protocol Handbook

Description

GeneAll® Allspin™ total DNA/RNA purification kit provides a convenient method for the isolation of total DNA and total RNA simultaneously from a single sample of tissue or cultured cells. DNA and RNA are purified separately from a same sample by individual but successive procedure using column B and column W respectively. Whole procedure can be performed in just 30 minutes and the length of obtained DNA is up to 50 kb (average is 30 kb) and that of RNA is longer than 200 nucleotides.

Features and Benefits

- Glassfiber membrane technology
- Sample size : Up to 20 mg tissue or up to 1×10^7 cultured cells
- Typical yield of RNA : Up to $20 \mu\text{g}$ / 1×10^6 cultured cells
Up to $60 \mu\text{g}$ / 10 mg liver tissue
- Typical yield of DNA : Up to $10 \mu\text{g}$ / 1×10^6 cultured cells
Up to $25 \mu\text{g}$ / 10 mg liver tissue
- High purity
- Preparation time : ~ 30 minutes
- Stable and consistent yield
- No phenol / chloroform extraction
- No ethanol precipitation
- Ready for use in RT-PCR, northern blotting, dot blotting, in vitro translation, molecular cloning, real-time PCR, RNase protection assays and other analytical procedures

The yield of genomic DNA and total RNA may vary depending on the tissue or cells from which it is obtained.

Material	Sample type	Average yield of genomic DNA	Average yield of Total RNA
Cultured cell ($\approx 1 \times 10^6$)	CHO	~ $7 \mu\text{g}$	~ $15 \mu\text{g}$
	RAW 264.7	~ $10 \mu\text{g}$	~ $20 \mu\text{g}$
Tissue (rat) (10 mg / prep)	Liver	~ $25 \mu\text{g}$	~ $60 \mu\text{g}$
	Kidney	~ $25 \mu\text{g}$	~ $30 \mu\text{g}$
	Brain	~ $12 \mu\text{g}$	~ $10 \mu\text{g}$
	Heart	~ $10 \mu\text{g}$	~ $9 \mu\text{g}$
	Spleen	~ $70 \mu\text{g}$	~ $80 \mu\text{g}$

For total RNA & DNA isolation from tissues and cultured cells

Allspin™
Column type B for DNA

Color : Red ring

Type : Spin

Maximum amount of starting samples : ~ 30 mg tissue or ~ 1×10^7 cellsMaximum loading volume : ~ 700 $\mu\ell$ Minimum elution volume : ~ 50 $\mu\ell$ Maximum binding capacity : ~ 100 μg

Nucleic acid binding size : ~ 50 kbp


Allspin™
Column type W for RNA

Blue ring

Spin

~ 30 mg tissue or ~ 1×10^7 cells~ 700 $\mu\ell$ ~ 30 $\mu\ell$ ~ 100 μg

> 200 nucleotides

Comparison Data


RT-PCR results from total RNA of rat heart tissue using Allspin™ and supplier kit were analysed on a 1% agarose gel.

M : 100 bp ladder

Lane 1, 2 : PCR of cDNA from Allspin™

Lane 3, 4 : PCR of cDNA from supplier A

Lane 5 : Negative control



PCR result from genomic DNA and total RNA eluate of CHO cells.

M : 100 bp ladder

Lane 1, 2 : Genomic DNA eluate from Allspin™


Lane 3, 4 : Genomic DNA eluate from Supplier A

Lane 5, 6 : Total RNA eluate from Allspin™

Lane 7, 8 : Total RNA eluate from Supplier A

Lane 9 : Negative control

Fig. I



Genomic DNA and total RNA were purified from RAW264.7 cells using Allspin™ and supplier A

Fig. I Genomic DNA were analysed on a 1% agarose gel

M : Lambda-HindIII

Lane 1, 2 : Genomic DNA from Allspin™

Lane 3, 4 : Genomic DNA from Supplier A

Fig. II.




Fig. II Total RNA were analysed on a 1% formaldehyde agarose gel

M : 0.5 ~ 10 kb RNA ladder

Lane 1, 2 : Total RNA from Allspin™

Lane 3, 4 : Total RNA from Supplier A

Products	Type	Size	Cat. No.
GeneAll® Allspin™	Spin	50	306-150

For RNA Purification

Hybrid-R™ / RiboEx™ / Ribospin™ / Allspin™ Series

RiboEx™ series are designed for total RNA isolation from various samples. RiboEx™ is based on the disruption of cells in a monophasic lysis solution containing phenol and salt followed by alcohol precipitation of the RNA. Hybrid-R™ eliminates alcohol precipitation by binding of RNA with column. RiboEx™ LS is a concentrated form of RiboEx™ and for total RNA isolation from liquid samples, while RiboEx™ is more suitable for solid samples and pelleted cells. Riboclear™ provides an easy and rapid method for RNA cleanup or concentration from various RNA samples in just 6 minutes. Ribospin™ series provide fast and easy method in convenient spin column format and isolate highly purified RNA in 15 minutes. Allspin™ total DNA / RNA purification kit provides a convenient method for the isolation of total DNA and total RNA simultaneously from a single sample of tissue or cultured cells.

	Hybrid-R™	Hybrid-R™ Blood RNA	Hybrid-R™ miRNA	RiboEx™	RiboEx™ LS	Ribospin™	Ribospin™ vRD (plus)*	Ribospin™ Plant	Riboclear™ (plus)***	Allspin™**
Specifications										
Format	Spin	Spin	Spin	Solution	Solution	Spin	Spin	Spin	Spin	Spin
Recommended sample volume	100 mg 1 x 10 ⁷ cells	250 µl	100 mg 1 x 10 ⁷ cells	100 mg 1 x 10 ⁷ cells	100 mg 250 µl	25 mg 1 x 10 ⁶ cells	300 µl	100 mg	100 µl	30 mg 1 x 10 ⁷ cells
Preparation time	30 min	30 min	30 min	50 ~ 65 min	50 ~ 65 min	15 min	20 min	25 min	6 min	30 min
Max. loading volume	700 µl	700 µl	700 µl	-	-	750 µl	800 µl	700 µl	800 µl	700 µl
Min. elution volume	30 µl	30 µl	30 µl	-	-	40 µl	30 µl	30 µl	30 µl	30 ~ 50 µl
Binding capacity	500 µg	100 µg	100 µg	-	-	100 µg	100 µg	100 µg	100 µg	100 µg

* Ribospin™ vRD *plus!* provides carrier RNA for purification of nucleic acid from very small amounts of sample.

** Allspin™ provides the method for the purification of genomic DNA and total RNA from tissues and cultured cells.

*** Riboclear™ *plus!* provides DNase for removal of DNA.

For RNA Purification

Hybrid-R™ / RiboEx™ / Ribospin™ / Allspin™ Series

	Hybrid-R™	Hybrid-R™ Blood RNA	Hybrid-R™ miRNA	RiboEx™	RiboEx™ LS	Ribospin™	Ribospin™ vRD (plus!)	Ribospin™ Plant	Riboclear™ (plus!)	Allspin™
Sample Type										
Animal cells	○	-	○	○	○	○	-	-	-	○
Animal tissues	○	-	○	○	△	○	-	-	-	○
Plant tissues	-	-	-	○	△	-	-	○	-	-
Bacteria	○	-	○	○	○	○	-	-	-	-
Yeast	○	-	○	○	○	-	-	-	-	-
Whole blood	-	○	-	-	○	-	-	-	-	-
Buffy coat	○	○	○	○	○	○	-	-	-	○
Various liquid sample	△	△	-	-	○	-	△	-	-	-
Viral sample	-	-	-	-	-	-	○	-	-	-
RNA cleanup/ concentration	-	-	-	-	-	-	-	-	○	-

○ Recommended / △ Suitable but not optimized



All columns in GeneAll® RNA related products are provided as individual packs (blister packs) to minimize the contamination.

GeneAll® AmpONE™

For more information about products, visit www.geneall.com

New! BB Solution!

BB solution

6X BB solution is newly developed gel loading buffer. And it's designed for easy loading and tracking of nucleic acids in agarose gels. In order to see PCR product or purified nucleic acids, it can be used as a loading dye. Even though PCR proceeds with BB solution, PCR works successfully and BB solution doesn't effect on the results. Furthermore PCR mixture can be loading directly without additional dye on agarose gels as PCR premix products.

gel loading dye

final conc. : 1X BB solution

example	
6X BB solution	1 μl
PCR product or DNA prep sample	5 μl

PCR mix

final conc. : 0.5X BB solution


example	
PCR reaction vol.	30 μl
6X BB solution	2.5 μl

PCR mixture (final conc. 0.5x)

final vol.	20 μl	30 μl	40 μl	50 μl	100 μl
6x BB	1.6 μl	2.5 μl	3.3 μl	4.1 μl	8.3 μl
PCR mix	18.4 μl	27.5 μl	36.7 μl	45.9 μl	91.7 μl

5. PCR Amplification System


AmpONE™ Taq DNA Polymerase	74
AmpONE™ α -Taq DNA Polymerase	75
AmpONE™ Pfu DNA Polymerase	76
AmpONE™ HS-Taq DNA Polymerase	77
AmpONE™ Clean Taq / Clean α -Taq DNA Polymerase	78
AmpONE™ Taq / α -Taq / HS-Taq Premix	79
AmpMaster™ Taq / α -Taq / HS-Taq	80
Selection Guide for DNA Polymerase	81



Comparison of PCR results using BB solution or house dye in RNA sample.
Sample : total RNA from CHO cells
Gel : 1% agarose gel

Sample A : 207 ng
Sample B : 253 ng

M : 1 kb ladder
Lane 1 : Add house dye to sample A
Lane 2 : Add house dye to sample B
Lane 3 : Add BB solution to sample A
Lane 4 : Add BB solution to sample B



Comparison of PCR results using BB solution or house dye for agarose gel loading
Template : human genomic DNA (50 ng)
Target size : 500 bp
Reaction vol. : 30 μ l
Loading vol. : 5 μ l

M1, M2 : 100 bp ladder
Lane 1, 2 : Using AmpONE™ Taq / α -Taq for PCR and house dye for gel loading
Lane 3, 4 : Using AmpONE™ Taq / α -Taq for PCR and using BB solution for gel loading
Lane 5, 6 : Using AmpONE™ Taq / α -Taq included BB solution and then directly gel loading

GeneAll® AmpONE™

Taq DNA Polymerase

AmpONE™ Taq DNA polymerase Reaction mix

for 50 μl reaction

10X Taq reaction buffer	5 μl
(optional : HQ buffer 5 ~ 20 μl)	
dNTP mix (2.5 mM each)	2 ~ 4 μl
primer 1	5 ~ 10 pmol
primer 2	5 ~ 10 pmol
template	0.1 ~ 100 ng
Taq (2.5 U/ μl)	0.5 ~ 1 μl
DW	up to 50 μl

Component list

Taq DNA Polymerase (2.5 U/ μl)
 10X Taq Reaction Buffer
 (with 25 mM Mg^{2+})
 dNTP Mix (2.5 mM each)
 HQ Buffer
 BB Solution
 Manual

Description


GeneAll® AmpONE™ Taq DNA polymerase is a recombinant enzyme derived from *Thermus aquaticus*, which is cloned and expressed in *E. coli* and possesses the same functions as the native enzyme. This enzyme is a thermostable DNA polymerase of 94 kDa and can be used in various experiments such as general PCR, RT-PCR and dideoxy-terminator-cycle sequencing. We have performed the quality control through activity test, purity test and endonuclease activity test.

Features and Benefits


- High fidelity, High purity
- Provides HQ buffer for the amplification of a higher order structure such GC-rich templates
- No 3' → 5' exonuclease activity
: addition of a single adenosine at 3' end of the extension product

Application


- General PCR
- TA-cloning
- DNA sequencing
- RT-PCR



Amplification of human genomic DNA.
 To check the amplification of various size the used primers are designed in various region.
 Template : human genomic DNA
 M : 1 kb ladder
 Lane 1 : 517 bp
 Lane 2 : 1.1 kb
 Lane 3 : 1.9 kb
 Lane 4 : 3.1 kb



Sensitivity of AmpONE™ Taq DNA polymerase on the quantity of template.
 Template : human genomic DNA
 M : 1 kb ladder
 Lane 1 : 20 ng
 Lane 2 : 10 ng
 Lane 3 : 1 ng
 Lane 4 : 100 pg
 Lane 5 : 10 pg



Comparison of AmpONE™ Taq DNA polymerase with other companies (1, 1.9 kb).
 Template : human gDNA (40 ng)
 M : 1 kb ladder
 Lane 1, 2 : GeneAll® AmpONE™ Taq
 Lane 3, 4 : company S
 Lane 5, 6 : company I
 Lane 7, 8 : company T

Products	Type	Size	Cat. No.
AmpONE™ Taq DNA Polymerase	(2.5U/ μl)	250 U	501-025
		500 U (250 U x 2)	501-050
		1000 U (250 U x 4)	501-100

dNTP mix and set are available. When ordering, please see the ordering information (page 82).

α-Taq DNA Polymerase

AmpONE™ α-Taq DNA polymerase Reaction mix

for 50 μℓ reaction

10X α-Taq reaction buffer	5 μℓ
(optional : HQ buffer 5 ~ 20 μℓ)	
dNTP mix (2.5 mM each)	2 ~ 4 μℓ
primer 1	5 ~ 10 pmol
primer 2	5 ~ 10 pmol
template	0.1 ~ 100 ng
α-Taq (2.5 U/μℓ)	0.5 ~ 1 μℓ
DW	up to 50 μℓ

Component list

α-Taq DNA Polymerase (2.5 U/μℓ)
 10X α-Taq Reaction Buffer
 (with 25 mM Mg²⁺)
 dNTP Mix (2.5 mM each)
 HQ Buffer
 BB Solution
 Manual

Description


GeneAll® AmpONE™ α-Taq DNA polymerase is a modified enzyme mixed Taq DNA polymerase with Pfu DNA polymerase which has proof-reading activity and has the ability to amplify a long PCR product (up to 20 kb). Although many other PCR enzymes with high fidelity, mainly derived from the *Pyrococcus furiosus* generally have a slow elongation rate, α-Taq DNA polymerase shows a fast elongation rate and more accurate PCR product formation.

Features and Benefits


- High fidelity, High purity
- Provides HQ buffer for the amplification of a higher order structure such GC-rich templates
- Addition of a single adenosine at 3' end of the extension product

Application


- General PCR
- Cloning for protein expression
- Long PCR not exceeding 20 Kb
- Multiplex PCR



Amplification of human genomic DNA.
 To check the amplification of various size the used primers are designed in various region.
 Template : human genomic DNA
 M : 1 kb ladder
 Lane 1 : 517 bp
 Lane 2 : 1.1 kb
 Lane 3 : 3.1 kb
 Lane 4 : 14 kb



Sensitivity of AmpONE™ α-Taq DNA polymerase on the quantity of template.
 Template : human genomic DNA
 M : 1 kb ladder
 Lane 1 : 20 ng
 Lane 2 : 10 ng
 Lane 3 : 1 ng
 Lane 4 : 100 pg
 Lane 5 : 10 pg



Comparison of AmpONE™ α-Taq DNA polymerase with other companies (3.8 kb).
 Template : human genomic DNA (40 ng)
 M : 1 kb ladder
 Lane 1, 2 : GeneAll® AmpONE™ α-Taq
 Lane 3, 4 : company S
 Lane 5, 6 : company I
 Lane 7, 8 : company T

Products	Type	Size	Cat. No.
AmpONE™ α-Taq DNA Polymerase	(2.5U/μℓ)	250 U	502-025
		500 U (250 U x 2)	502-050
		1000 U (250 U x 4)	502-100

dNTP mix and set are available. When ordering, please see the ordering information (page 82).

Pfu DNA Polymerase

AmpONE™ Pfu DNA polymerase Reaction mix

for 50 µℓ reaction

- 10X Pfu reaction buffer 5 µℓ
(optional : HQ buffer 5 ~ 20 µℓ)
- dNTP mix (2.5 mM each) 2 ~ 4 µℓ
- primer 1 5 ~ 10 pmol
- primer 2 5 ~ 10 pmol
- template 1 ~ 100 ng
- Pfu (2.5 U/µℓ) 0.5 ~ 1 µℓ
- DW up to 50 µℓ

Component list

- Pfu DNA Polymerase (2.5 U/µℓ)
- 10X Pfu Reaction Buffer (with 25 mM Mg²⁺)
- dNTP Mix (2.5 mM each)
- HQ Buffer
- BB Solution
- Manual

Description



AmpONE™ Pfu DNA polymerase is a recombinant modified enzyme derived from Pyrococcus furiosus, which is cloned and expressed in E. coli. This enzyme has 3' → 5' exonuclease (proofreading) activity as well as polymerase activity, providing higher fidelity than Taq polymerase. And it is especially suited to PCR for protein expression and site direct mutagenesis.

Features and Benefits


- High fidelity, High purity
- Provides HQ buffer for the amplification of a higher order structure such GC-rich templates

Application


- Cloning for protein expression
- Site direct mutagenesis
- Blunt-end cloning



Amplification of human genomic DNA.
To check the amplification of various size the used primers are designed in various region.
Template : human genomic DNA
M : 1 kb ladder
Lane 1 : 517 bp
Lane 2 : 1.1 kb
Lane 3 : 1.9 kb
Lane 4 : 3.1 kb



Sensitivity of AmpONE™ Pfu DNA polymerase on the quantity of template.
Template : human genomic DNA
M : 1 kb ladder
Lane 1 : 20 ng
Lane 2 : 10 ng
Lane 3 : 1 ng
Lane 4 : 100 pg
Lane 5 : 10 pg



Comparison of AmpONE™ Pfu DNA polymerase with other companies.
Template : human genomic DNA (40 ng)
M : 1 kb ladder
Lane 1, 2 : GeneAll® AmpONE™ Pfu
Lane 3, 4 : company S (cloned type)
Lane 5, 6 : company S (native type)
Lane 7, 8 : company I

Products	Type	Size	Cat. No.
AmpONE™ Pfu DNA Polymerase	(2.5U/µℓ)	250 U	503-025
		500 U (250 U x 2)	503-050
		1000 U (250 U x 4)	503-100

dNTP mix and set are available. When ordering, please see the ordering information (page 82).

HS-Taq DNA Polymerase

AmpONE™ HS-Taq DNA polymerase Reaction mix

for 50 μ l reaction

10X HS-Taq reaction buffer	5 μ l
(optional : HQ buffer 5 ~ 20 μ l)	
dNTP mix (2.5 mM each)	4 μ l
primer 1	5 ~ 10 pmol
primer 2	5 ~ 10 pmol
template	0.1 ~ 100 ng
HS Taq (2.5 U/ μ l)	0.5 ~ 1 μ l
DW	up to 50 μ l

Component list

HS-Taq DNA Polymerase (2.5 U/ μ l)
 10X HS-Taq reaction buffer
 (with 25 mM Mg²⁺)
 dNTP mix (2.5 mM each)
 HQ Buffer
 BB Solution
 Manual

Description


AmpONE™ HS-Taq DNA Polymerase, a modified form of **GeneAll® AmpONE™ Taq DNA Polymerase**, is designed to enhance the specificity, sensitivity and yield of DNA amplification. The activity of AmpONE™ HS-Taq DNA Polymerase is rapidly restored during the initial denaturation step of PCR. By limiting polymerase activity prior to PCR cycling, the amplification of non-specifically annealed primers or primer dimers is reduced and target yield is increased.

Features and Benefits

- Reduced nonspecific amplification
- High specificity
- Enhanced sensitivity
- Convenient PCR set-up at room temperature
- Provides HQ buffer for the amplification of a higher order structure such GC-rich templates

Application

- Highly specific amplification
- Real-time PCR
- RT-PCR
- Multiplex PCR
- TA-cloning



Products	Type	Size	Cat. No.
AmpONE™ HS-Taq DNA Polymerase	(2.5U/ μ l)	250 U	531-025
		500 U (250 U x 2)	531-050
		1000 U (250 U x 4)	531-100

dNTP mix and set are available. When ordering, please see the ordering information (page 82).

GeneAll® AmpONE™

Clean Taq / Clean α-Taq DNA Polymerase

AmpONE™ Clean Taq DNA polymerase Reaction mix

for 50 μℓ reaction

- 10X Clean Taq reaction buffer 5 μℓ
(optional : HQ buffer 5 ~ 20 μℓ)
- dNTP mix (2.5 mM each) 4 μℓ
- primer 1 5 ~ 10 pmol
- primer 2 5 ~ 10 pmol
- template 0.1 ~ 100 ng
- Clean Taq (2.5 U/μℓ) 0.5 ~ 1 μℓ
- DW up to 50 μℓ

AmpONE™ Clean α-Taq DNA polymerase Reaction mix

for 50 μℓ reaction

- 10X Clean α-Taq reaction buffer 5 μℓ
(optional : HQ buffer 5 ~ 20 μℓ)
- dNTP mix (2.5 mM each) 4 μℓ
- primer 1 5 ~ 10 pmol
- primer 2 5 ~ 10 pmol
- template 0.1 ~ 100 ng
- Clean α-Taq (2.5 U/μℓ) 0.5 ~ 1 μℓ
- DW up to 50 μℓ

Description


AmpONE™ Clean Taq and Clean α-Taq DNA Polymerase are thermostable enzymes to minimize the contamination of genomic DNA using new purification method. There is no E. coli genomic DNA contamination in the amplification (over 40 cycles) using E. coli 16S RNA primer.

Features and Benefits

- Genomic DNA contamination FREE
- Enhanced sensitivity
- Providing HQ buffer (PCR enhancer buffer) for the amplification of a higher order structure
- Applications : Long PCR using AmpONE™ Clean α-Taq (~ 20 Kb)
real-time PCR, RT-PCR, Multiplex PCR, TA-cloning

Component list

- 10X Clean Taq / Clean α-Taq DNA Polymerase reaction buffer (with 25 mM Mg²⁺)
- dNTP Mix (2.5 mM each)
- HQ Buffer
- BB Solution
- Manual



Genomic DNA contamination test using E. coli 16S RNA primer
M : Marker, G : GeneAll® AmpONE™ Clean Taq, S : company S
T : company T, I : company I

Products	Type	Size	Cat. No.
AmpONE™ Clean Taq DNA Polymerase	(2.5U/μℓ)	250 U	551-025
		500 U (250 U x 2)	551-050
		1000 U (250 U x 4)	551-100
AmpONE™ Clean α-Taq DNA Polymerase	(2.5U/μℓ)	250 U	552-025
		500 U (250 U x 2)	552-050
		1000 U (250 U x 4)	552-100

dNTP mix and set are available. When ordering, please see the ordering information(page 82).

Taq / α -Taq / HS-Taq PremixAmpONE[™] Premix
Reaction mix

<Lyophilized>

Reaction vol.	20 μ l	50 μ l
Premix	~ μ l	~ μ l
primer 1 (10 pmole / μ l)	1 μ l	1 ~ 2 μ l
primer 2 (10 pmole / μ l)	1 μ l	1 ~ 2 μ l
template	1 ~ 50 ng	1 ~ 100 ng
DW	up to 20 μ l	up to 50 μ l
final reaction vol.	20 μ l	50 μ l

<Solution>

Reaction vol.	20 μ l	50 μ l
Premix	10 μ l	25 μ l
primer 1 (10 pmole / μ l)	1 μ l	1 ~ 2 μ l
primer 2 (10 pmole / μ l)	1 μ l	1 ~ 2 μ l
template	1 ~ 50 ng	1 ~ 100 ng
DW	up to 20 μ l	up to 50 μ l
final reaction vol.	20 μ l	50 μ l

Component list


Taq / α -Taq / HS-Taq premix 96 tubes
(0.2 ml 8-tube strip x 12ea)
0.2 ml PCR tube storage rack
Manual

Description

AmpONE[™] Taq / α -Taq / HS-Taq Premix contains all reaction components required for PCR, such as reaction buffer, dNTP, gel loading dye, stabilizer and sediment in addition to Taq DNA polymerase. It is recommended to use in routine PCR (below 10 kb), TA cloning and primer extension. This mixture is stable for 1 year at -20 °C or 2 weeks at room temperature. It is ready-to-use mixture pipetting steps are minimized, reducing the possibility of errors and contamination. Room temperature reaction setup using this mixture is fast and easy. Included loading dye migrates through 1.0% agarose gels run in 0.5X TBE at approximately the same rate as DNA 300 bp in length.

Features and Benefits

- Lyophilized mix type or 2X solution type
- Use of low tracking dye (1% agarose gel, 300 bp)
- Ready to use
- Minimal handling
- Stable for 1 year at -20 °C or 2 weeks at RT
- Offer of PCR tube rack



Amplification of AmpONE[™] Taq Premix.
To check the amplification of various size the used primers are designed in various region.
Template : human genomic DNA (27 ng)
M : 1 kb ladder
Lane 1 : 514 bp
Lane 2 : 1 kb
Lane 3 : 1.9 kb
Lane 4 : 3.8 kb
Reaction vol. : 20 μ l Loading vol. : 2 μ l

Products	Format	Type	Size	Cat. No.
AmpONE [™] Taq Premix	Lyophilized	20 μ l	96 tubes	521-200
		50 μ l	96 tubes	521-500
AmpONE [™] α -Taq Premix	Lyophilized	20 μ l	96 tubes	522-200
		50 μ l	96 tubes	522-500
AmpONE [™] Taq Premix (w/o dye)	Lyophilized	20 μ l	96 tubes	524-200
AmpONE [™] HS-Taq Premix	Solution	20 μ l	96 tubes	525-200
		50 μ l	96 tubes	525-500
AmpONE [™] Taq Premix	Solution	20 μ l	96 tubes	526-200
		50 μ l	96 tubes	526-500
AmpONE [™] α -Taq Premix	Solution	20 μ l	96 tubes	527-200
		50 μ l	96 tubes	527-500

GeneAll® AmpMaster™

Taq / α-Taq / HS-Taq

AmpMaster™ Reaction mix

Reaction vol.	20 μℓ	50 μℓ
2X Master mix	10 μℓ	25 μℓ
primer 1 (10 pmole / μℓ)	1 μℓ	1 ~ 2 μℓ
primer 2 (10 pmole / μℓ)	1 μℓ	1 ~ 2 μℓ
template	1 ~ 50 ng	1 ~ 100 ng
DW	up to 20 μℓ	up to 50 μℓ
final reaction vol.	20 μℓ	50 μℓ

Component list

2X Taq / α-Taq / HS-Taq PCR master mix
Manual


Description

GeneAll® AmpMaster™ series contain all reaction components required for PCR, such as reaction buffer, dNTP, gel loading dye, stabilizer and sediment in addition to Taq / α-Taq / HS-Taq DNA polymerase. It is recommended for use in routine PCR (below 10 kb), TA cloning and primer extension.

GeneAll® AmpMaster™ are stable for 1 year at -20 °C or 2 months at 4 °C. It is ready-to-use mixture pipetting steps are minimized, reducing the possibility of errors and contamination. Room temperature reaction setup using this mixture is fast and easy. Included loading dye migrates through 1.0% agarose gels run in 0.5X TBE at approximately the same rate as DNA 300 bp in length.

Features and Benefits

- 2X solution type
- Use of low tracking dye (1% agarose gel, 300 bp)
- Easy reaction setup
- Fewer pipetting steps
- Stable for 1 year at -20 °C or 2 months at 4 °C



Consistency test of AmpMaster™ Taq.
Template : human genomic DNA (27 ng)
M : 1 kb plus ladder
Lane 1 ~ 8 : 1 kb
Reaction vol. : 20 μℓ Loading vol. : 2 μℓ

Products	Type	Size	Cat. No.
AmpMaster™ Taq	1 ml	0.5 ml x 2 tubes	541-010
	5 ml	0.5 ml x 10 tubes	541-050
AmpMaster™ Taq (w/o dye)	1 ml	0.5 ml x 2 tubes	544-010
	5 ml	0.5 ml x 10 tubes	544-050
AmpMaster™ α-Taq	1 ml	0.5 ml x 2 tubes	542-010
	5 ml	0.5 ml x 10 tubes	542-050
AmpMaster™ HS-Taq	1 ml	0.5 ml x 2 tubes	545-010
	5 ml	0.5 ml x 10 tubes	545-050

For PCR Amplification System

AmpONE™ / AmpMaster™ Series

GeneAll® AmpONE™ DNA polymerase series are designed for various applications such as cloning, genotyping, sequencing, routine PCR, real-time PCR and long PCR and provided with HQ buffer for the amplification of a higher order structure. AmpONE™ Premix (lyophilized or solution form) and AmpMaster™ mix (liquid form) are provided with DNA polymerase, reaction buffer, dNTP and loading dye. The BB solution contains sediment and two dyes (xylene cyanol FF and Orange G) that will allow its location in the gel to be identified and separated individual strands to be monitored during electrophoresis. This solution can be applied to use of two purposes. The first case is to use as a general gel loading dye. And the second is for master mix in pre-PCR. In this case, PCR products can be loaded directly into the wells of the gel without additional loading dye.

	AmpONE™ Taq DNA Polymerase	AmpONE™ α -Taq DNA Polymerase	AmpONE™ Pfu DNA Polymerase	AmpONE™ HS-Taq DNA Polymerase
Applications				
Standard PCR (< 3 kb)	+++	+++	+	++
Hot start PCR	-	-	-	+++
Multiplex PCR	+	+++	+	++
Nested PCR	+	++	+	+
Amplification product size	< 5 kb	< 20 kb	< 5 kb	< 5 kb
HQ Buffer*	○	○	○	○
BB Solution**	○	○	○	○


• **GeneAll® AmpONE™** Premix and Master mix are available.

• **GeneAll® AmpONE™** Premix and Master mix are made from AmpONE™ Taq / α -Taq / HS-Taq DNA polymerase which contain all reaction components required for PCR, such as reaction buffer, dNTP, gel loading dye, stabilizer and sediment.

• **GeneAll® AmpONE™** Premix is lyophilized or solution form and **GeneAll® AmpMaster™** mix is 2-fold concentrated form.

* A novel additive that enables efficient amplification of GC-rich template or long size amplification.

** High-quality gel-loading dye for analysis of PCR samples using electrophoresis.




Comparison of Taq DNA polymerase with other companies (1, 1.9 kb)

Template : human genomic DNA (40 ng)

M : 1 kb ladder


Lane 3, 4 : company S

Lane 7, 8 : company T



Lane 1, 2 : **GeneAll® AmpONE™** Taq

Lane 5, 6 : company I




Comparison of α -Taq DNA polymerase with other companies (14, 20 kb)

Template : human genomic DNA (100 ng)

M : Lambda-HindIII

Lane 3, 4 : company S

Lane 7, 8 : company T



Lane 1, 2 : **GeneAll® AmpONE™** α -Taq

Lane 5, 6 : company I

GeneAll[®] SV Column & Filter Column Accessories (Mini)



Column type Q



Column type B



Column type Qe



Column type W



Column type D



Column type F



Column type G



Column type V



Column type S



EzClear[™] filter Column



Micro Column type S



EzSep[™] filter Column

Ordering Information

Products	Scale	Size	Cat. No.	Type
GeneAll® Hybrid-Q™ for rapid preparation of plasmid DNA				
Plasmid Rapidprep		50	100-150	mini / spin
		100	100-102	
GeneAll® Exprep™ for purification of plasmid DNA				
Plasmid SV	mini	50	101-150	spin / vacuum
		200	101-102	
		1,000	101-111	
	Midi	26	101-226	spin / vacuum
		50	101-250	
		100	101-201	
GeneAll® Exfection™ for purification of highly pure plasmid DNA				
Plasmid LE (Low Endotoxin)	mini	50	111-150	spin / vacuum
		200	111-102	
	Midi	26	111-226	spin / vacuum
100		111-201		
Plasmid EF (Endotoxin Free)	Midi	20	121-220	spin
		100	121-201	
GeneAll® Expin™ for purification of fragment DNA				
Gel SV	mini	50	102-150	spin / vacuum
		200	102-102	
PCR SV	mini	50	103-150	spin / vacuum
		200	103-102	
CleanUp SV	mini	50	113-150	spin / vacuum
		200	113-102	
Combo GP	mini	50	112-150	spin / vacuum
		200	112-102	
GeneAll® Exgene™ for isolation of total DNA				
Tissue SV	mini	100	104-101	spin / vacuum
		250	104-152	
	Midi	26	104-226	spin / vacuum
		100	104-201	
	MAXI	10	104-310	spin / vacuum
		26	104-326	
Tissue <i>plus!</i> SV	mini	100	109-101	spin / vacuum
		250	109-152	
	Midi	26	109-226	spin / vacuum
		100	109-201	
	MAXI	10	109-310	spin / vacuum
		26	109-326	

Ordering Information

Products	Scale	Size	Cat. No.	Type
GeneAll® Exgene™ for isolation of total DNA				
Blood SV	mini	100	105-101	spin / vacuum
		250	105-152	
	Midi	26	105-226	spin / vacuum
		100	105-201	
	MAXI	10	105-310	spin / vacuum
		26	105-326	
Cell SV	mini	100	106-101	spin / vacuum
		250	106-152	
	MAXI	10	106-310	spin / vacuum
		26	106-326	
Clinic SV	mini	100	108-101	spin / vacuum
		250	108-152	
	Midi	26	108-226	spin / vacuum
		100	108-201	
	MAXI	10	108-310	spin / vacuum
		26	108-326	
Genomic DNA micro		50	118-050	spin
Plant SV	mini	100	117-101	spin / vacuum
		250	117-152	
	Midi	26	117-226	spin / vacuum
		100	117-201	
	MAXI	10	117-310	spin / vacuum
		26	117-326	
Soil DNA mini	mini	50	114-150	spin
GMO SV	mini	50	107-150	spin / vacuum
		200	107-102	
GeneAll® GenEx™ for isolation of total DNA				
GenEx™ Blood	Sx	100	220-101	solution
		500	220-105	
GenEx™ Cell	Lx	100	220-301	solution
	Sx	100	221-101	solution
		500	221-105	
GenEx™ Tissue	Lx	100	221-301	solution
		100	221-301	
	GenEx™ Plant	Sx	100	222-101
500			222-105	
Lx		100	222-301	solution
GenEx™ Plant	Sx	100	227-101	solution
		Mx	100	
	Lx	100	227-301	
GenEx™ Plant <i>plus!</i>	Sx	100	228-101	solution
		Mx	50	
	Lx	20	228-320	

※ For more information about ordering, visit www.geneall.com

Products	Scale	Size	Cat. No.	Type
----------	-------	------	----------	------

GeneAll® DirEx™ Single tube DNA preparation buffer for PCR

DirEx™		50	250-050	solution
--------	--	----	---------	----------

GeneAll® RNA series for purification of total RNA

RiboEx™	mini	100	301-001	solution
		200	301-002	
Hybrid-R™	mini	100	305-101	spin
Hybrid-R™ Blood RNA	mini	50	315-150	spin
Hybrid-R™ miRNA	mini	50	325-150	spin
RiboEx™ LS	mini	100	302-001	solution
		200	302-002	
Riboclear™	mini	50	303-150	spin
Riboclear™ <i>plus!</i>	mini	50	313-150	spin
Ribospin™	mini	50	304-150	spin
Ribospin™ vRD	mini	50	302-150	spin
Ribospin™ vRD <i>plus!</i>	mini	50	312-150	spin
Ribospin™ Plant	mini	50	307-150	spin
Allspin™	mini	50	306-150	spin

GeneAll® AmpONE™ for PCR amplification

Taq DNA polymerase		250 U	501-025	(2.5 U/μl)
		500 U	501-050	
		1,000 U	501-100	
α- Taq DNA polymerase		250 U	502-025	(2.5 U/μl)
		500 U	502-050	
		1,000 U	502-100	
Pfu DNA polymerase		250 U	503-025	(2.5 U/μl)
		500 U	503-050	
		1,000 U	503-100	
Hotstart Taq DNA polymerase		250 U	531-025	(2.5 U/μl)
		500 U	531-050	
		1,000 U	531-100	
Clean Taq DNA polymerase		250 U	551-025	(2.5 U/μl)
		500 U	551-050	
		1,000 U	551-100	
Clean α- Taq DNA polymerase		250 U	552-025	(2.5 U/μl)
		500 U	552-050	
		1,000 U	552-100	

Ordering Information

Products	Scale	Size	Cat. No.	Type
GeneAll® AmpONE™ for PCR amplification				
Taq Premix	96 tubes	20 µl	521-200	lyophilized
		50 µl	521-500	
		20 µl	526-200	solution
		50 µl	526-500	
α - Taq Premix	96 tubes	20 µl	522-200	lyophilized
		50 µl	522-500	
		20 µl	527-200	solution
		50 µl	527-500	
HS - Taq Premix	96 tubes	20 µl	525-200	solution
		50 µl	525-500	
Taq Premix (w/o dye)	96 tubes	20 µl	524-200	lyophilized
dNTP Mix		500 µl	509-020	2.5 mM each
dNTP set (set of dATP, dCTP, dGTP and dTTP)		1 ml x 4 tubes	509-040	100 mM
GeneAll® AmpMaster™ for PCR amplification				
Taq Master mix		2x	541-010	0.5 ml x 2 tubes
		2x	541-050	0.5 ml x 10 tubes
α - Taq Master mix		2x	542-010	0.5 ml x 2 tubes
		2x	542-050	0.5 ml x 10 tubes
HS - Taq Master mix		2x	545-010	0.5 ml x 2 tubes
		2x	545-050	0.5 ml x 10 tubes

* Each dNTP is available

Visit GeneAll® Community

www.geneall.com
www.geneall.co.kr

Mailing Address

GeneAll Biotechnology Co., LTD.
GeneAll Bldg., 303-7 Dongnam-ro,
Songpa-gu, Seoul, Korea 138-859

Ordering information

Tel: 82-2-407-0096
Fax: 82-2-407-0779
E-mail: sales@geneall.com

Technical information

Tel: 82-2-407-0096
Fax: 82-2-407-0779
E-mail: tech@geneall.com

Customer & Technical Support

Do not hesitate to ask us any question.
We thank you for any comment or advice.

