

2X HyperScript™ One-step RT-PCR master mix

Ver. 1.00

Cat. No. 602-110 (0.5ml x 2 tubes)

Storage at : -20 °C

Lot. No. M0202RPM

Expiration date : 2014.07.30

Description

Hyperscript™ one-step RT-PCR master mix is a 2x premix ready to use for Reverse Transcriptase(RT) reaction and Polymerase Chain Reaction(PCR). This mix contains Hyperscript RTase and AmpONE™ HS-TaQ DNA polymerase and both RT and PCR reactions are carried out successively in a single tube. HS-TaQ DNA polymerase remains inactivated until RT reaction is completed, and it is turned on at high temperature of PCR cycle. HS Taq polymerase can amplify the fragment up to 5 kb in length. The reaction volume can be adjusted according to the experimental purpose.

Hyperscript™ one-step RT-PCR master mix contains all reaction components required for RT and PCR, such as reaction buffer, dNTPs, RNase inhibitor and stabilizer in addition to enzymes, except primers and templates.

Components

2X HyperScript™ one-step RT-PCR mix tubes (0.5 ml x 2ea)

Storage condition

Stable for 1 year at -20°C

Ingredients of Hyperscript™ one-step RT-PCR master mix

(based on 20 µl reaction volume)

Hyperscript™ Reverse Transcriptase

HS-TaQ DNA polymerase

dNTPs (mixture)

Reaction buffer (2.5 mM MgCl₂)

Stabilizer

RNase inhibitor

Procedure

- RNA template

The following 20-µl reaction volume can be used for 1 pg – 2 µg of total RNA or 10 pg – 500 ng of poly(A) RNA.

- One step RT-PCR

1. Combine the followings in nuclease-free tube

Components	Concentration	Amount
One-step RTPCR master mix	2X	10 µl
Total RNA		x µl (up to 2 µg)
Forward Primer	10 pmol/µl	1 µl
Reverse Primer	10 pmol/µl	1.5 µl
DEPC treated water		Adjust the mixture volume to 20 µl with water

2. Add each component, mix gently and collect by brief centrifugation. Do reaction as follows.

Step	Temp.	Time	Cycle No.
cDNA synthesis step	42 - 60°C (recommend- 55°C)	30 - 60 min	1 cycle
Pre-Denaturation	94°C	2 - 5 min	1 cycle
Denaturation	94°C	30 - 60 sec	25- 35 cycles
Annealing	X°C	30 - 60 sec	
Extension	72°C	1 min / kb	
Post-Extension	72°C	2-5 min	1 cycle

* The reaction temperature can be increased to the extent of 65°C if the target gene has a somewhat intricate secondary structure.