

HyperScript™ One-step RT-PCR premix

Ver. 1.00

Cat. No. 602-102

Storage at : -20 °C

Lot. No. M0201RPP

Expiration date : 2014.05.30

Description

HyperScript™ one-step RT-PCR premix is a 2x premix ready to use for Reverse Transcriptase(RT) reaction and Polymerase Chain Reaction(PCR). The component of this premix is optimized for 20 $\mu\ell$ of reaction volume which is widely used for many downstream applications. This mix contains HyperScript™ RTase and AmpONE™ HS-Taq DNA polymerase and both RT and PCR reactions are carried out successively in a single tube. HS-Taq DNA polymerase remains inactivated until RT reaction is completed, and it is turned on at high temperature of PCR cycle. HS-Taq polymerase can amplify the fragment up to 5 kb in length.

HyperScript™ one-step RT-PCR premix contains all reaction components required for RT and PCR, such as reaction buffer, dNTPs, RNase inhibitor and stabilizer in addition to enzymes, except primers and templates.

Components

HyperScript™ one-step RT-PCR premix 96 tubes (8-tubes(0.2 ml) strip x 12ea with PCR tube rack)

Storage condition

Stable for 1 year at -20°C

Ingredients of HyperScript™ one-step RT-PCR premix

(based on 20 $\mu\ell$ reaction volume)

HyperScript™ Reverse Transcriptase

HS-Taq DNA polymerase

dNTPs (mixture)

Reaction buffer (2.5 mM MgCl₂)

Stabilizer

RNase inhibitor

Procedure

(1) RNA template

The following 20- $\mu\ell$ reaction volume can be used for 1 pg – 2 μg of total RNA or 10 pg – 500 ng of poly(A) RNA.

(2) One step RT-PCR

1. Add the followings into 8-strip tubes

Components	Concentration	Amount
One step RT-PCR premix	2X	10 $\mu\ell$
Add the followings		
Total RNA		x $\mu\ell$ (up to 2 μg)
Forward Primer	10 pmol/ $\mu\ell$	1 $\mu\ell$
Reverse Primer	10 pmol/ $\mu\ell$	1.5 $\mu\ell$
DEPC treated water		Adjust the mixture volume to 20 $\mu\ell$ with water

2. Add each component, mix gently and collect by brief centrifugation. Perform reaction as follows.

Step	Temp.	Time	Cycle No.
cDNA synthesis step	42 - 60°C (recommend- 55°C)	30 - 60 min	1 cycle
Pre-Denaturation	94°C	2 - 5 min	1 cycle
Denaturation	94°C	30 - 60 sec	25- 35 cycles
Annealing	X °C	30 - 60 sec	
Extension	72°C	1 min / kb	
Post-Extension	72°C	2-5 min	1 cycle

* The reaction temperature can be increased to the extent of 65°C if the target gene has a somewhat intricate secondary structure.