

# Ribospin™ vRD

VIRAL RNA/DNA PURIFICATION HANDBOOK

REF

302-150/302-103



HB3200



GeneAll Bldg., 303-7 Dongnam-ro Songpa-gu,  
Seoul, South Korea 138-859

GeneAll

CE  
IVD

## **Customer & Technical Support**

Should you have any further questions, do not hesitate to contact us.  
We appreciate your comments and advice.

## **Contact Information**

[www.geneall.com](http://www.geneall.com)

Tel : 82-2-407-0096

Fax : 82-2-407-0779

E-mail(Order/Sales) : [sales@geneall.com](mailto:sales@geneall.com)

E-mail(Tech. Info.) : [tech@geneall.com](mailto:tech@geneall.com)

## **Visit GeneAll® Community**

[www.geneall.com](http://www.geneall.com)

[www.geneall.co.kr](http://www.geneall.co.kr)

## **Trademarks**

Hybrid-Q™, Exprep™, Exfection™, Expin™, Exgene™, GenEx™, DirEx™, Hybrid-R™, RiboEx™, Ribospin™,  
Riboclear™, Allspin™, RiboSaver™, EzClear™, EzSep™, EzPure™, EzPass™, AmpONE™, AmpMaster™,  
RealAmp™, HyperScript™, ProtinEx™, PAGESTA™, STEADI™ are trademarks of GeneAll Biotechnology Co.,  
Ltd.

© 2017 GeneAll Biotechnology, all rights reserved.

This protocol handbook is included in :

GeneAll® Ribospin™ vRD (302-150, 203-103)

Visit [www.geneall.com](http://www.geneall.com) or [www.geneall.co.kr](http://www.geneall.co.kr) for FAQ, Q&A and more information.

# Used symbols and Markings

<b>REF</b>	Catalog number	<b>IVD</b>	In vitro diagnostic medical device
<b>LOT</b>	Batch number	<b>H B</b>	Handbook code
	Use by		Consult instruction for use
	Manufacturer information		Contains sufficient for <N> tests
	Do not reuse		Temperature limitation
	Production date	<b>EC REP</b>	European Authorized Representative
	Important note	<b>CONC</b>	Contains the concentrated solution. Additional material must be added before use
	Write down the current date after adding ethanol to the bottle	<b>EtOH ?</b> <input checked="" type="checkbox"/>	Mark up after adding ethanol

# INDEX

- 04**      Index
- 05**      Kit Contents
- Product Specifications
- 06**      Quality Control
- Storage Conditions
- Safety Information
- Preventing RNase Contamination
- 07**      Product Description
- 08**      Protocol
- 10**      Troubleshooting Guide
- 12**      Ordering Information

## Kit Contents

Components	Quantity		Storage
	Cat. No.	302-150	302-103
No. of preparation	50	300	
Buffer VL	30 ml	170 ml	
Buffer RB I (concentrate) *	8 ml	48 ml	
Buffer RBW (concentrate) *	13 ml	77 ml	
Buffer RNW (concentrate) * †	6 ml	34 ml	
Nuclease-free water	15 ml	20 ml	Room temperature (15~25°C)
Mini column type V (with collection tube)	50	300	
1.5 ml microcentrifuge tube	50	300	
Protocol Handbook	1	1	



\* Before first use, add absolute ethanol (ACS grade or better) into buffer RB I, RBW and RNW as indicated on the bottle.

† Contains sodium azide as a preservative

## Product Specifications

### Ribospin™ vRD

Type	Using spin column
Maximum volume of starting samples	300 µl/prep
Preparation time	20 min
Maximum loading volume	800 µl
Minimum elution volume	30 µl

## Quality Control

All components of GeneAll® Ribospin™ vRD are manufactured in strictly clean conditions, and their degree of cleanliness is monitored periodically.

To maintain consistency, a quality control process is carried out thoroughly from lot to lot and only the qualified kits are approved for delivery according to ISO 9001:2008 and EN ISO13485:2012.

## Storage Conditions

All components of GeneAll® Ribospin™ vRD should be stored at room temperature (15~25°C).

During shipment or storage under cool ambient condition, a precipitate can form in buffer VL. In such a case, heat the bottle to 56°C to dissolve completely. GeneAll® Ribospin™ vRD is guaranteed until the expiration date printed on the product box.

## Safety Information

The buffers included in the GeneAll® Ribospin™ vRD contain irritants which are harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken when handling such materials. Always wear gloves and eye protection, and follow standard safety precautions.



Buffer VL, RB I, and RBW contain chaotropes, which can form highly reactive compounds when combined with bleach. Do NOT add bleach or acidic solutions directly to the sample-preparation waste.

## Preventing RNase contamination

RNase can be introduced accidentally during RNA purification. Wear disposable gloves always, because skin often contains bacteria and molds that can be a source of RNase contamination. Use sterile, disposable plastic wares and automatic pipettes to prevent cross-contamination of RNase from shared equipment.

## Product Description

GeneAll® Ribospin™ vRD provides a convenient method for isolation of RNA and DNA from cell-free fluid, cell-culture medium, plasma, serum, swab, urine and virus-infected samples.

GeneAll® Ribospin™ vRD utilizes the glass fiber membrane technology for the fastest and the most convenient nucleic acid isolation as a sufficient level for downstream application instead of conventional alcohol precipitation or phenol/chloroform extraction.

The buffer system of Ribospin™ vRD provides the effective binding condition of RNA and DNA to glass fiber membrane and the impurities on the membrane are washed away by two different wash buffers. At least, pure RNA and DNA are eluted in Nuclease-free water. The whole procedure may take only 15 minutes at room temperature and the eluate is suitable for PCR, RT-PCR, or any downstream application without further manipulation. The purified nucleic acid should be treated with care because RNA is very sensitive to contaminants such as RNases, often found on general lab ware and dust. To ensure RNA-stability after extraction, it is recommended to store at 4°C for immediate analysis or to freeze at -70°C for long-term storage.

# PROTOCOL FOR

## Ribospin™ vRD

### Equipment and reagents to be supplied by user

- \* Ethanol (>99%, ACS grade or better)
- \* 1.5 ml microcentrifuge tubes
- \* Micropipettes and sterile pipet tips
- \* Centrifuge capable of attaining 10,000 x g
- \* Vortex mixer



- Ethanol (>99%, ACS grade or better) must be added before the first use of buffer RB I, RBW and RNW. Please refer to the information on the label of each bottle.
- If a precipitate is formed in buffer VL, heat to 56°C to dissolve completely before use.

- 1. Transfer up to 300 µl sample (swab-storage media, cell-free fluid, cell-culture supernatant, plasma, serum, urine) in 1.5 ml microcentrifuge tube.**
- 2. Add 500 µl buffer VL to the tube and lyse the sample by pipetting or vortexing.**

The volume of buffer VL can be adjusted in proportion to the volume of sample.  
For proper lysis, the complete mixing of sample and buffer VL is essential.
- 3. Incubate the lysate for 10 min at room temperature.**

After this step, briefly centrifuge the tube to remove drops from the inside of the lid.
- 4. Add 700 µl buffer RB I to the lysate and mix thoroughly by inverting or vortexing.**

The volume of buffer RB I can be adjusted in proportion to the volume of lysate.  
\* Do NOT centrifuge at this step.
- 5. Transfer up to 750 µl of the mixture to a mini column type V.**

**6. Centrifuge at  $\geq 10,000 \times g$  for 30 sec at room temperature.**

Discard the pass-through and reinsert the mini column back into the same tube.

**7. Repeat step 5~6 with the remainder of the sample.**

Discard the pass-through and reinsert the mini column back into the same tube.

**8. Add 500  $\mu l$  buffer RBW to the mini column.**

**9. Centrifuge at  $\geq 10,000 \times g$  for 30 sec at room temperature.**

Discard the pass-through and reinsert the mini column back into the same tube.

**| 10. Add 500  $\mu l$  buffer RNW to the mini column.**

**| 11. Centrifuge at  $\geq 10,000 \times g$  for 30 sec at room temperature.**

Discard the pass-through and reinsert the mini column back into the same tube.

**| 12. Centrifuge at  $\geq 10,000 \times g$  for an additional 1 min at room temperature to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube (provided).**

Residual ethanol may interfere with downstream reactions. Care must be taken at this step for eliminating the carry-over buffer RNW.

If the carry-over buffer RNW still occurs, centrifuge again for 1 min at full speed before transferring the column to the new 1.5 ml microcentrifuge tube.

**| 13. Add 30~50  $\mu l$  of Nuclease-free water to the center of the membrane in the mini column.**

**Let it stand for 1 min.**

**| 14. Centrifuge at  $\geq 10,000 \times g$  for 1 min at room temperature.**

Purified nucleic acid can be stored at 4°C for immediate analysis or at -70°C for long-term storage.

## Troubleshooting Guide

Facts	Possible Causes	Suggestions
<b>Low yield</b>	<b>Poor quality of starting material</b>	Fresh sample or well-conserved sample should be used for good result. Repeated freezing and thawing the sample should be avoided.
	<b>Low concentration of viral particle in the starting sample</b>	Use more starting sample. If the amount of sample is more than 300 $\mu$ l, concentrate the volume to 300 $\mu$ l using a micro-concentrator.
	<b>Inefficient or insufficient lysis</b>	Be sure to incubate for 10 minutes at room temperature after adding buffer VL. For proper lysis, the complete mixing of the sample and buffer VL is essential.
	<b>Improper elution</b>	Add Nuclease-free water to the center of the mini column membrane and perform incubation for 1 minute before centrifugation.
	<b>Precipitate in buffer VL</b>	A precipitate can be formed in buffer VL at cool ambient temperature. It is because the buffer VL is saturated and its solubility would be reduced at low temperature. Before experiment, any precipitate in the buffer VL should be dissolved completely by heating the buffer at 56°C or above until it disappears.
	<b>Degradation of RNA</b>	RNase can be introduced during purification of nucleic acid. Be certain not to introduce any RNases during the procedure of later handling. Keep tubes closed whenever possible during the extraction and use RNase-free products with sterile and disposable plastic ware.
	<b>Buffer RBI, RBW, or RNW was prepared incorrectly</b>	Check that the concentrated buffer RBI, RBW, and RNW were diluted with the correct volume of absolute ethanol.

## Troubleshooting Guide

Facts	Possible Causes	Suggestions
<b>Purified nucleic acid does not perform well in down-stream application</b>	<b>Residual ethanol from buffer RNW remains in eluate</b>	Care must be taken for eliminating the carry-over buffer RNW before elution step. The membrane of mini column should be kept completely dry via additional centrifugation (Step 12, page 9) or air-drying.
	<b>Incorrect order of buffer RBW and RNW</b>	Ensure that buffer RBW and RNW are used in the correct order during extraction. If used in the wrong order, perform the last washing step with buffer RNW.

# Ordering Information

Products	Scale	Size	Cat. No.	Type	Products	Scale	Size	Cat. No.	Type
<b>GeneAll® Hybrid-Q™ for rapid preparation of plasmid DNA</b>									
Plasmid Rapidprep	mini	50	100-150						
		200	100-102	spin					
<b>GeneAll® Exprep™ for preparation of plasmid DNA</b>									
	mini	50	101-150	spin /					
		200	101-102	vacuum					
Plasmid SV	Midi	26	101-226						
	Midi	50	101-250	spin /					
		100	101-201	vacuum					
<b>GeneAll® Exfection™ for preparation of transfection-grade plasmid DNA</b>									
	mini	50	111-150	spin /					
Plasmid LE (Low Endotoxin)		200	111-102	vacuum					
	Midi	26	111-226	spin /					
	Midi	100	111-201	vacuum					
Plasmid EF (Endotoxin Free)	Midi	20	121-220						
		100	121-201	spin					
<b>GeneAll® Expi™ for purification of fragment DNA</b>									
Gel SV	mini	50	102-150	spin /					
		200	102-102	vacuum					
PCR SV	mini	50	103-150	spin /					
		200	103-102	vacuum					
CleanUp SV	mini	50	113-150	spin /					
		200	113-102	vacuum					
Combo GP	mini	50	112-150	spin /					
		200	112-102	vacuum					
<b>GeneAll® Exgene™ for isolation of total DNA</b>									
	mini	100	104-101	spin /					
		250	104-152	vacuum					
Tissue SV	Midi	26	104-226	spin /					
		100	104-201	vacuum					
	MAXI	10	104-310	spin /					
		26	104-326	vacuum					
	mini	100	109-101	spin /					
		250	109-152	vacuum					
Tissue plus! SV	Midi	26	109-226	spin /					
		100	109-201	vacuum					
	MAXI	10	109-310	spin /					
		26	109-326	vacuum					
<b>GeneAll® Exgene™ for isolation of total DNA</b>									
	mini	100	105-101	spin /					
		250	105-152	vacuum					
Blood SV	Midi	26	105-226	spin /					
		100	105-201	vacuum					
	MAXI	10	105-310	spin /					
		26	105-326	vacuum					
Cell SV	mini	100	106-101	spin /					
		250	106-152	vacuum					
	MAXI	10	106-310	spin /					
		26	106-326	vacuum					
Clinic SV	mini	100	108-101	spin /					
		250	108-152	vacuum					
	Midi	26	108-226	spin /					
		100	108-201	vacuum					
	MAXI	10	108-310	spin /					
		26	108-326	vacuum					
Genomic DNA micro		50	118-050	spin					
	mini	100	117-101	spin /					
		250	117-152	vacuum					
Plant SV	Midi	26	117-226	spin /					
		100	117-201	vacuum					
	MAXI	10	117-310	spin /					
		26	117-326	vacuum					
Soil DNA mini	mini	50	114-150	spin					
Stool DNA mini	mini	50	115-150	spin					
Viral DNA / RNA	mini	50	128-150	spin					
FFPE Tissue DNA	mini	50	138-150	spin					
		250	138-152	spin					
<b>GeneAll® Exgene™ for isolation of total DNA</b>									
	Sx	100	220-101						
		500	220-105						
GenEx™ Blood	Lx	100	220-301						
	Sx	100	221-101						
		500	221-105						
GenEx™ Cell	Lx	100	221-301						
	Sx	100	222-101						
		500	222-105						
GenEx™ Tissue	Lx	100	222-301						

Products	Scale	Size	Cat. No.	Type	Products	Scale	Size	Cat. No.	Type	
<b>GeneAll® GenEx™ for isolation of total DNA</b>										
GenEx™ Plant	Sx	100	227-101	solution	Taq DNA polymerase	250 U	501-025			
	Mx	100	227-201			500 U	501-050	(2.5 U/ $\mu$ l)		
	Lx	100	227-301			1,000 U	501-100			
GenEx™ Plant plus!	Sx	100	228-101	solution	$\alpha$ -Taq DNA polymerase	250 U	502-025			
	Mx	50	228-250			500 U	502-050	(2.5 U/ $\mu$ l)		
	Lx	20	228-320			1,000 U	502-100			
<b>GeneAll® DirEx™ series</b> for preparation of PCR-template without extraction										
DirEx™		100	250-101	solution	$\alpha$ -Pfu DNA polymerase	250 U	504-025			
DirEx™ Fast-Tissue		96 T	260-011	solution		500 U	504-050	(2.5 U/ $\mu$ l)		
DirEx™ Fast-Cultured cell		96 T	260-021	solution		1,000 U	504-100			
DirEx™ Fast-Whole blood		96 T	260-031	solution	Fast-Pfu DNA polymerase	250 U	505-025			
DirEx™ Fast-Blood stain		96 T	260-041	solution		500 U	505-050	(2.5 U/ $\mu$ l)		
DirEx™ Fast-Hair		96 T	260-051	solution		1,000 U	505-100			
DirEx™ Fast-Buccal swab		96 T	260-061	solution	Hotstart Taq DNA polymerase	250 U	531-025			
DirEx™ Fast-Cigarette		96 T	260-071	solution		500 U	531-050	(2.5 U/ $\mu$ l)		
						1,000 U	531-100			
<b>GeneAll® RNA series</b> for preparation of total RNA										
RiboEx™	mini	100	301-001	solution	Taq Premix	20 $\mu$ l	521-200	lyophilized		
		200	301-002			50 $\mu$ l	521-500			
Hybrid-R™	mini	100	305-101	spin		20 $\mu$ l	526-200	solution		
Hybrid-R™ Blood RNA mini		50	315-150	spin		50 $\mu$ l	526-500			
Hybrid-R™ miRNA	mini	50	325-150	spin	$\alpha$ -Taq Premix	20 $\mu$ l	522-200	lyophilized		
RiboEx™ LS	mini	100	302-001	solution		50 $\mu$ l	522-500			
		200	302-002			20 $\mu$ l	527-200			
Riboclear™	mini	50	303-150	spin		50 $\mu$ l	527-500			
Riboclear™ plus!	mini	50	313-150	spin	HS-Taq Premix	20 $\mu$ l	525-200	solution		
Ribospin™	mini	50	304-150	spin		50 $\mu$ l	525-500			
Ribospin™ II	mini	50	314-150	spin		20 $\mu$ l	520-200	lyophilized		
		300	314-103			50 $\mu$ l	520-500			
Ribospin™ vRD	mini	50	302-150	spin	$\alpha$ -Pfu Premix	20 $\mu$ l	523-500	solution		
Ribospin™ vRD plus!	mini	50	312-150	spin		50 $\mu$ l	524-200	lyophilized		
Ribospin™ vRD II	mini	50	322-150	spin	dNTPs mix	500 $\mu$ l	509-020	2.5 mM each		
Ribospin™ Plant	mini	50	307-150	spin	dNTPs set	1 ml x 4 tubes	509-040	100 mM		
Ribospin™ Seed / Fruit	mini	50	317-150	spin						
Allspin™	mini	50	306-150	spin						
RiboSaver™	mini	100	351-001	solution						

Products	Scale	Size	Cat. No.	Type	Products	Size	Cat. No.
<b>GeneAll® AmpMaster™ for PCR amplification</b>							
Taq Master mix	0.5 ml x 2 tubes	541-010	solution		ProteinEx™	100 ml	701-001
	0.5 ml x 10 tubes	541-050	solution		Animal cell / tissue		solution
α-Taq Master mix	0.5 ml x 2 tubes	542-010	solution		PAGESTA™		
	0.5 ml x 10 tubes	542-050	solution		Reducing	1 ml x 10 tubes	751-001
HS-Taq Master mix	0.5 ml x 2 tubes	545-010	solution		5X SDS-PAGE		solution
	0.5 ml x 10 tubes	545-050	solution		Sample Buffer		
α-Pfu Master mix	0.5 ml x 2 tubes	543-010	solution				
	0.5 ml x 10 tubes	543-050	solution				
<b>GeneAll® HyperScript™ for Reverse Transcription</b>							
Reverse Transcriptase	10,000 U	601-100	solution		<b>GeneAll® STEADI™ for automatic nucleic acid purification</b>		
RT Master mix	0.5 ml x 2 tubes	601-710	solution		STEADI™ 12 Instrument		GST012
RT Master mix with oligo (dT) <sub>20</sub>	0.5 ml x 2 tubes	601-730	solution		STEADI™ 24 Instrument		GST024
RT Master mix with random hexamer	0.5 ml x 2 tubes	601-740	solution		STEADI™ Genomic DNA Cell / Tissue	96	401-104
RT Premix	96 tubes, 20 µl	601-602	solution		STEADI™ Genomic DNA Blood	96	402-105
RT Premix with oligo (dT) <sub>20</sub>	96 tubes, 20 µl	601-632	solution		STEADI™ Bacteria DNA	96	403-106
RT Premix with random hexamer	96 tubes, 20 µl	601-642	solution		STEADI™ Total RNA	96	404-304
One-step RT-PCR Master mix	0.5 ml x 2 tubes	602-110	solution		STEADI™ Viral DNA / RNA	96	405-322
One-step RT-PCR Premix	96 tubes, 20 µl	602-102	solution		STEADI™ CFC Seed DNA / RNA	96	406-C02
First strand Synthesis Kit	50 reaction	605-005	solution				
ZymAll™ RNase Inhibitor	10,000 U	605-010	solution				
ZymAll™ RNase Inhibitor	4,000 U	605-004	solution				
<b>GeneAll® RealAmp™ for qPCR amplification</b>							
SYBR qPCR Master mix (2X, Low ROX)	200 rxn 20 µl	801-020	solution				
	500 rxn 20 µl	801-050					
SYBR qPCR Master mix (2X, High ROX)	200 rxn 20 µl	801-021	solution				
	500 rxn 20 µl	801-051					

**NOTE**

---



**GENEALL BIOTECHNOLOGY CO., LTD**

GeneAll Bldg., 303-7 Dongnam-ro,  
Songpa-gu, Seoul, South Korea 138-859  
E-mail : sales@geneall.com  
[www.geneall.com](http://www.geneall.com)

Tel : 82-2-407-0096

Fax : 82-2-407-0779

©2017 GeneAll Biotechnology, All rights reserved

EC REP

MT Promedt Consulting GmbH  
Altenhofstrasse 80, 66386 St. Ingbert,  
Germany  
TEL : +49 (0) 6894 581020  
FAX : +49 (0) 6894 581021

Edited by **JBR**  
Designed by **KDY**  
2017.02